everStick™ PERIO

Fibre reinforcement for periodontal splinting

DISTRIBUTED BY

GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium
TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A., TEL: +1-708-597-0900
www.gc.dental/america

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529541, TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia, TEL: +61 2 9301 8200

GC SOUTH AMERICA
Rua Heliodora, 399, Santana - São Paulo, SP, Brasil, TEL.:+55-11-2925-0965
CEP: 02022-051, CNPJ: 08.279.999/0001-61, RESP. TÊC: Erick de Lima - CRO/SP 100.866

RESPONSIBLE MANUFACTURER IN CANADA

GC AMERICA INC.
3737 W. 127th Street, Alsip, IL 60803 U.S.A.

Fibre type: Silanated E-Glass fibre impregnated with Bis-GMA and PMMA

Form: Unidirectional fibre bundle

Diameter: ~ 1.1-1.3 mm

Keep away
from sunlight

Rx Only

UK
CA
0086

GC EUROPE AG
ZÜRICHSTRASSE 31,
6004 LUZERN SWITZERLAND

Caution: US Federal law restricts this device to sale by or on the order of a dentist.

30000532-NU6813


MADE IN FINLAND

everStick™ PERIO

- EN FIBRE REINFORCEMENT FOR PERIODONTAL SPLINTING
- DE GLASFASERVERSTÄRKUNG FÜR PARODONTALE SCHIENUNGEN
- FR FIBRE DE RENFORT POUR CONTENTION PARODONTALE
- IT RINFORZO IN FIBRA PER SPLINTAGGIO PARODONTALE
- ES FIBRA DE REFUERZO PARA FERULIZACIONES PERIODONTALES
- NL VEZELVERSTERKING VOOR PERIRADICULAIR SPALKEN
- DA FIBERFORSTÆRKET SPLINTING
- SV FIBERFÖRSTÄRKT SPLINTING
- PT FIBRA PARA REFORÇO DE CONTENÇÕES PERIODONTAIS
- EL ΙΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΓΙΑ ΕΝΙΣΧΥΣΗ ΠΕΡΙΟΔΟΝΤΙΚΗΣ ΑΚΙΝΗΤΟΠΟΙΗΣΗΣ
- FI LASIKUITULUJITE KISKOTUKSIIN
- NO FIBERFORSTERKET SPLINTING
- ZH 复合树脂纤维用于牙周夹板治疗
- ID PENGUAT SERAT UNTUK SPLINTING PERIODONTAL

EN	Fibre reinforcement for periodontal splinting	5	SV	Fiberförstärkt splinting	51
DE	Glasfaserverstärkung für parodontale Schienungen	11	PT	Fibra para reforço de contenções periodontais	57
FR	Fibre de renfort pour contention parodontale	18	EL	ΪΝΕΣ ΥΑΛΟΝΗΜΑΤΩΝ ΓΙΑ ΕΝΙΣΧΥΣΗ ΠΕΡΙΟΔΟΝΤΙΚΗΣ ΑΚΙΝΗΤΟΠΟΙΗΣΗΣ	64
IT	Rinforzo in fibra per splintaggio parodontale	24	FI	Lasikuitulujite kiskotuksiin	71
ES	Fibra de refuerzo para ferulizaciones periodontales	31	NO	Fiberforsterket splinting	77
NL	Vezelversterking voor periradicaalr spalken	38	ZH	用于牙周夹板治疗的复合树脂纤维	84
DA	Fiberforstærket splinting	46	ID	Penguat Serat untuk Splinting Periodontal	89

For use only by a dental professional in the indications for use.

PACKAGES

everStickPERIO 2 x 12 cm

everStickPERIO 1 x 8 cm

everStickINTRO

8 cm everStickC&B; 8 cm everStickPERIO;
1x StickREFIX D silicone instrument

everStickCOMBI

8cm everStickC&B; 8cm everStickPERIO;
30 cm² everStickNET; 5x everStickPOST Ø1.2;
1x StickSTEPPER hand instrument;
1x StickREFIX D silicone instrument

everStickSTARTER KIT

8 cm everStickC&B; 8 cm everStickPERIO;
30cm² everStickNET; 5 x everStickPOST
Ø 1.2; 6 ml GC Modeling Liquid bottle;
2 ml G-ænial Universal Flo syringe A2;
20 dispensing tips and 1 light protective cap;

1x StickREFIX D
silicone instrument; 1x StickSTEPPER;
1x StickCARRIER

WHAT IS everStickPERIO?

everStickPERIO fibre reinforcement is a combination of glass fibres and a permeable polymer/resin gel matrix for use in dentistry as a reinforcing material. The polymer/resin gel holds the individual glass fibres in a bundle, which facilitates handling of the fibres. The fibre bundle is flexible and sticky, which allows it easily and reliably bond to teeth. The primary application for everStickPERIO fibre reinforcement is the splinting of teeth.

INDICATIONS FOR USE

everStickPERIO fibre reinforcement is recommended for use as periodontal surface-retained and intracoronal splints

- Lingual/palatal splints
- Labial splints

- Occlusal splints

CONTRAINDICATION

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

COMPOSITION

Silane treated e-type glass fiber roving, methacrylates, initiators, inhibitors.

INSTRUCTIONS FOR USE:

IMPORTANT: The everStickPERIO fibres should be positioned close to the incisal edge to minimise the forces the splint will be subjected to. Also, the splint should not interfere with the occlusal contacts e.g., on the palatal surface of the upper anterior teeth.

Spot fixation with composite does not provide a sufficient bond between the fibre

splint and tooth surface. Bond the fibre splint to the teeth for their entire length. Cover the fibre bundle with a thin (0.5 mm) layer of composite, including the approximal areas, when bonding it to the tooth surfaces. At the occlusal contact in intra-coronal splints, the optimal thickness of the composite layer on top of the fibre is approximately 1–2 mm.

Splints in the upper anterior area do not necessarily have any room palatally for a surface-retained splint, due to the occlusion. Consequently, an intra-coronal grooved splint or surface retained labial splint must be considered.

In an anterior area's lingual/palatal splint, a short additional fibre can be used to offset the occlusal forces that tend to loosen fibre splints in the canine region. The additional fibre is attached to the labial surface of the canine and the lateral incisor.

SURFACE-RETAINED PERIODONTAL SPLINTING IN ANTERIOR AREA

1. Measure and cut the fibre

Measure the length of fibre needed from the dental arch by using, for example, a periodontal probe or dental floss, to prepare an everStickPERIO fibre splint. Open the foil package and use tweezers to pull out an appropriate amount of silicone embedded fibres. Using sharp scissors cut the appropriate amount of fibre needed along with the silicone. Shield the fibre from light by placing it under a cover during preparation of the teeth to be bonded. Close the foil bag tightly with its sticker. Keep the bag in a refrigerator (at a temperature of 2-8°C/35.6-46.4°F) when you do not need it.

2. Clean the tooth surfaces

The entire length of the fibre splint must be bonded to tooth surfaces. Clean the tooth surfaces with a paste of pumice and water, rinse and air-dry the area. Place wedges in

the approximal spaces as necessary, so that the spaces to be cleaned are not filled with composite. If you are working without wedges, be careful not to block these spaces with composite – see item 5.

3. Etch the tooth surfaces

Etch the tooth surfaces and inter-proximal spaces thoroughly with ortho-phosphoric acid, in the area of the splint, in accordance with the instructions of the bonding agent manufacturer. Preferably etch slightly wider than necessary rather than too little. The recommended enamel etching time for surface-retained areas is 45 to 60 seconds. Rinse with water and air-dry the tooth surfaces thoroughly after etching. As with all bonded restorations a dry operating field is necessary and rubber dam isolation is highly recommended.

4. Bond the tooth surfaces

Use the adhesive bonding technique for bonding teeth according to the instructions of the bonding agent manufacturer.

Apply the bonding agent to the entire area to be bonded. Light-cure the bonding agent as described by the manufacturer.

5. Apply flowable composite

Apply a thin layer of flowable composite (for example, G-ænial Universal Flo or G-ænial Universal Injectable) on to the surface of the teeth for the entire width of the fibre bundle. Carefully cover the bonding area with a thin layer (about 0.5 mm) of composite including the approximal spaces. Leave enough space for cleaning the approximal spaces. Do not cure the composite during this phase.

6. Position and light-cure the fibre

Remove the white protective paper and use tweezers to pick the fibre up from the silicone groove. Remove any residual silicone from the fibre bundle. Position the fibre bundle on top of the uncured flowable composite. Aim to place the fibre as incisally as possible in the anterior area. Make sure that it will not be in occlusion.

Position one end of the fibre bundle first by pressing it down with the StickSTEPPER or other handinstrument. Pre-cure the fibre in place, one tooth at a time, for about five seconds, using a curing light.

The wide-tipped StickSTEPPER or other handinstrument shields the rest of the fibre from light. As the fibre is made of light-conducting material, it is recommended to direct the light-curer away from the uncured fibre bundle. Press the fibre into the approximal spaces as well. Make sure that the embrasures are not blocked with fibre and composite.

7. Cover and finish the splint

After pre-curing, cover the entire fibre splint with a thin layer of composite. Note that in surface-retained areas the fibre bundle can be coated with a thin layer (0.5 mm) of composite. Then light cure the whole splint for 40 seconds, one tooth or coverage area of the light-curer at a time. Be careful not to cut the fibre when finishing/polishing the

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The steps for the intra-coronal splint are the same as for the surface retained splint except for the preparation of a groove, mesial to distal, in the teeth to be splinted.

1. Prepare the groove for the teeth

Prepare a groove for the teeth to be splinted with a minimum width of 2 mm. Ideally, it should stay within the enamel as this will provide the best bond. At the occlusal contact, the optimal thickness of the composite layer on top of the fibre is 1–2 mm. The margins of the groove should be bevelled as this will ensure the best marginal integrity and enlarge the etched enamel surface for composite bonding. Place wedges in the approximal spaces, if possible, so that the embrasures are not filled with composite. If you are working without wedges, note that

the embrasures must remain free of composite.

2. Measure and cut the fibre

Measure the length of fibre needed to prepare an everStickPERIO fibre splint for the prepared groove by using, for example, a periodontal probe or dental floss. Open the foil package and use tweezers to pull out an appropriate amount of silicone embedded fibres. Using sharp scissors cut the appropriate amount of fibre needed along with the silicone. Shield the fibre from light by placing it under a cover during preparation of the teeth to be bonded. Close the foil bag tightly with its sticker. Keep the bag in a refrigerator (at a temperature of 2-8°C/35.6-46.4F) when you do not need it.

3. Etch the teeth

Etch the prepared groove thoroughly with ortho-phosphoric acid in accordance with the instructions of the bonding agent manufacturer. Rinse with water and air-dry the tooth surfaces thoroughly after etching.

As with all bonded restorations a dry operating field is necessary and rubber dam isolation is highly recommended.

4. Bond the teeth

Bond the prepared groove area in accordance with the instructions of the bonding agent manufacturer. Light-cure the bonding agent as described by the manufacturer.

5. Apply flowable composite

Apply a thin layer of flowable composite (for example, G-ænial Universal Flo or G-ænial Universal Injectable) into the prepared groove. Be careful not to block the embrasures with composite. Do not light cure the composite during this phase.

6. Position and light cure the fibre

Remove the white protective paper and use tweezers to pick the fibre up from the silicone groove. Remove any residual silicone from the fibre bundle. Position the fibre bundle on top of the uncured flowable composite. Aim to place the fibre as incisally as possible in the anterior area.

Make sure that it will not be in occlusion. Position one end of the fibre bundle first by pressing it down with the StickSTEPPER or other hand instrument. Pre-cure the fibre in place, one tooth at a time, for about five seconds, using a curing light. The wide-tipped StickSTEPPER or other hand instrument shields the rest of the fibre from light. As the fibre is made of light-conducting material, it is recommended to direct the light-curer away from the uncured fibre bundle. Press the fibre into the approximal spaces as well. Make sure that the embrasures are not blocked with fibre and composite.

7. Cover and finish the splint

Cover the entire fibre splint and fill the groove with a 1-2 mm composite layer. Then light cure the whole splint for 40 seconds, one tooth or coverage area of the light-curer at a time. Be careful not to cut the fibre when finishing/polishing the splint.

TIPS AND RECOMMENDATIONS

- The fibre must always be covered entirely with composite.
- Use rubber dam isolation to ensure a dry operating field.
- Always try to handle the fibres with instruments to avoid contamination e.g. powdered gloves.
- Use the StickSTEPPER or other handinstrument for sectional curing of the fibre bundle.
- Position the splint as close as possible to the incisal edges of the teeth in the anterior area.
- If, after placing the fibre, you notice it is too long, shorten it with a diamond bur during the finishing phase of the splint. Apply a light curable unfilled and solvent-free methacrylate resin (e.g. GC Modeling Liquid) to the exposed fibre surface, blow the resin in to a thin layer and light cure. Cover the fibre again carefully with composite.
- Avoid cutting the fibres during the finishing phase.

- At the occlusal contact area, the optimal thickness of the composite layer on top of the fibre splint is approximately 1–2 mm.

STORING: everStick products should always be stored in a refrigerator (2-8°C/35.6-46.4°F). In addition, the products should be protected from light by packing them in the sealed foil package after use. An elevated temperature and exposure to bright light may shorten the lifetime of everStick products.

Prior to application, the products are taken out of the refrigerator and the foil package opened, but kept away from bright daylight or artificial light. While cutting the fibre bundle, the rest of the fibre bundle inside the foil package should be kept covered from light. Immediately after cutting a sufficient length for the fibre construction, the foil package is carefully resealed and returned to the refrigerator.

TIPS AND RECOMMENDATIONS

1. The everStick products should be used clinically with care and the patient should be warned not to abrade the fitting surface to avoid exposing irritation-causing fibres.
2. If the surface of the fibre bundle feels dry, but it is fully bendable and not polymerized, adding a drop of a light curable unfilled and solvent-free methacrylate resin (e.g., GC Modeling Liquid) will return the flexibility/workability of the material. Polymerization can be observed as white spots at bending area when bending the bundle.
3. The everStick fibres do not achieve their full strength immediately after the final light-curing of 40 seconds. The polymerization of the fibres will continue during the next 24 hours.

CAUTION

4. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. The use of powder

free gloves is recommended with everStick products.

5. Unpolymerised resin can cause skin sensitisation to acrylates in some people. If your skin comes in contact with an encounter resin, wash it thoroughly with soap and water. Avoid contact of uncured material with skin, mucous membrane, or eyes. Unpolymerised everStick products may have a slight irritating effect and lead to sensitization to methacrylates in rare cases.

6. Polymerize everStick before waste disposal.

7. Do not use the product if the primary package of aluminium foil pouch is damaged. Product may be pre-polymerized and not usable.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
<https://www.gc.dental/europe>.

They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental

UNDESIRE EFFECTS- REPORTING:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system:
vigilance@gc.dental

In this way you will contribute to improve the

safety of this product.

UK Responsible Person
GC UNITED KINGDOM Ltd.
Coopers Court Newport Pagnell
Buckinghamshire
MK16 8JS
United Kingdom

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Revisi terakhir : 05/2024

EN Undesired effects - Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

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