GC Fuji TRIAGE™ CAPSULE

RADIOOPAKE GLASS IONOMER PROTECTION AND TEMPORARY RESTORATIVE MATERIAL IN CAPSULES

For use only by a dental professional in the indications for use.

- INDICATIONS FOR USE**
- Fluoride protection.
 - Root surface protection.
 - Hyperplastic prevention and control.
 - Protection for mature enamel.
 - Temporary filling of endodontic access.
 - Endosseous restorations.

- CONTRAINDICATIONS**
- Direct pulp capping.
 - Avoid use of this product in patients with known glass ionomer cement allergies.

COMPOSITION
Powder: Fluoro-alumino-silikatna pasta, pigment
Liquid: Water, polyalkylacrylate, carboxylic acid

DIRECTIONS FOR USE
PINK Shade is constant. WHITE Shade is chemical set only

Powder / Liquid Ratio (g/g)	0.30 / 0.15
Mixing Time (sec)	140*
Working Time (23°C/73°F) (min. sec.)	140*
Working Time (37°C/99°F) (min. sec.)	730*
Finishing Time (37°C/99°F, from start of mixing) (min. sec.)	600*
Finishing Time / Light cure (Pink Command set 37°C/99°F, from start of mixing) (min. sec.)	400*

ISO 9917-1: Glass Protection Cement, Restoration

APPLICATION

ROOT SURFACE PROTECTION

HYPERPLASTIC PREVENTION AND CONTROL PROTECTION FOR IMMATURE ENAMEL

- Preparation of the tooth surfaces (e.g. for fluoride protection or root surface protection).
- After cleaning the tooth surfaces (prophylaxis with pumice and water) in usual manner, clean the tooth surface with water. Avoid aggraving the overcured enamel.
- Apply the capsule directly on the tooth surface (prophylactic and sealant in water) to obliterate any gaps, promote adhesion of water. Remove water by the vacuum suction.
- The conditioner should be rinsed off thoroughly after application.
- By dry blotting with a cotton pellet or gently blowing with an air syringe (Fig. A-1). DO NOT DESICCATE. Avoid moisture in contact with the tooth surface.
- Apply the capsule directly on the tooth surface. Avoid immediately seal it into a mixer (or an amalgamator) for use for 10 seconds (Fig. A-2).

- Placement: Immediately remove the mixed capsule to the mixer and lead into the GC CAPSULE APPLIER. Do not allow to remain in the capsule. Dispense (Fig. A-3). The working time is 1 minute 40 seconds from the start of mixing at 23°C (73°F). Higher temperature will shorten the working time.
- Extrude the mixture onto the tooth surface (Fig. A-4) then a brush or spread a thin film of GC Fuji TRIAGE CAPSULE directly over the root surface of hyperplastic area or over the occlusal surface.
- Use a faster set is desired, use a visible light curing device for 20-40 seconds. Place light source as closely as possible to the tooth surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.
- To adjust the direction of the nozzle, hold the nozzle with the capsule towards you and rotate the capsule body.
- To remove the used capsule, push the applicator release button. Tilt the capsule body upwards.
- After placement, when the material starts to lose the adhesive appearance (or after curing with the light curing device), apply GC Fuji VARNISH (blow dry or GC Fuji COAT LC light cure) to the sealed area and the margins using a cotton pellet or sponge.
- After curing under air/water spray can be performed 6 minutes from start mix (chemically set) or 4 minutes from start mix (photo-curing) for a silicone finishing mix (Fig. A-8).
- Apply GC Fuji VARNISH or GC Fuji COAT LC (Fig. A-9).

TEMPORARY FILLING OF ENDODONTIC ACCESS

- Cleaning the pulp chamber.
- After adequate pulp chamber clean and gently dry the pulp chamber with an air syringe (Fig. B-1).
- Fill the chamber with a cotton pellet (Fig. B-2).
- See directions in A (above), Section 2.
- See directions in A (above), Section 3, a and b.
- Extrude the mixture directly over the cotton pellet (Fig. B-3).
- Use a faster set is desired, use a visible light curing device* for 20-40 seconds. Place light source as closely as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.
- Moisture protection and Finishing. See directions in A, Section 3, a and f.

INTERMEDIATE RESTORATIONS

- Remove casts and loose debris with hand instruments.
- Apply a thin film of gently clean the carious surface with CAVITY CONDITIONER for 10 seconds or DENTIN CONDITIONER for 20 seconds.
- For better retention of the cement, dry by blotting with a cotton pellet or gently blowing with an air syringe (Fig. C-1). DO NOT DESICCATE. Avoid moisture in contact with the tooth surface.
- Apply the capsule directly on the tooth surface (Fig. C-2). Form a uniform film and cover the tooth.
- Placement: Immediately remove the mixed capsule to the mixer and lead into the GC CAPSULE APPLIER. Do not allow to remain in the capsule. Dispense (Fig. C-3). The working time is 1 minute 40 seconds from the start of mixing at 23°C (73°F). Higher temperature will shorten the working time.
- Extrude the mixture directly over the prepared lesion or tooth surface (Fig. C-4). Then a brush or spread a thin film of GC Fuji TRIAGE CAPSULE directly over the lesion or tooth surface.
- Use a faster set is desired, use a visible light curing device* for 20-40 seconds. Place light source as closely as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.
- Moisture protection and Finishing. See directions in A, Section 3, a and f.

SHADE

Pink, White

STORAGE

Recommended for optimal performance, store at temperature of 2-25°C (36-77°F).

PACKAGE

Average contents per capsule: 0.30 g powder and 0.15 g (0.12 ml) liquid. Minimum net volume of material per capsule: 0.13 ml.

CAUTION

- In case of contact with oral tissue or skin, rinse immediately with a copious amount of water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In case cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gccorp.com>

They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUADMED data sheets (<https://ec.europa.eu/euadmed/>) or contact us at regulatory@gc.com

Undesired effects: Reporting. If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this IFU, please report them directly through the relevant contact person, by selecting the proper authority of your country accessible through http://www.gc.europa.eu/wh/contacts/medical-devices-contacts_en as well as to our internal vigilance system: vigilance@gc.com

This way you will contribute to improve the safety of this product.

Last revised: 03/2024

MANUFACTURED BY GC CORPORATION
1-1-1 Higashi, Isehara, Saitama, Tokyo 147-8685, Japan

DISTRIBUTED BY GC CORPORATION
1-1-1 Higashi, Isehara, Saitama, Tokyo 147-8685, Japan

EU: GC EUROPE NV
Pleinlaan 29, 1830 Louvain-la-Neuve, Belgium
Tel: +32 (0) 21 40 10 00

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