30000236 - NU5105
IFU G-ænial SHADE GUIDE

G-ænial™ SHADE GUIDE

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

Used for shade selection.

PRODUCT DESCRIPTION

G-ænial is a light-cure, radio-opaque restorative material to be used intra-orally and classified as a Type 1 and Class 2 (Group 1) per ISO standard 4049.

G-ænial ANTERIOR has a radiopacity equivalent to 1,5 - 2,0 mm of aluminum (dentine = 1 mm, enamel = 2 mm) and , G-ænial POSTERIOR has a radiopacity equivalent to 2,5 - 3,0 mm of aluminum and TE shade has no radiopacity.

The particle size of inorganic fillers range is 0,02 - 2,9 µm for the ANTERIOR and TE shade, 0,02 - 1,0 µm for the POSTERIOR.

The total amount of inorganic filler is approximately 63 vol% for the ANTERIOR, 59 vol% for TE shade, 65 vol% for the POSTERIOR.

HOW TO USE

1. Clean the tooth with pumice and water. Shade selection should be made prior to isolation. Keep the teeth hydrated. Select the appropriate G-ænial shades by referring to the G-ænial shade guide.

2. Each shade tab consists of four levels with 1,0, 1,5, 2,0 and 3,0 mm in thickness respectively. When layered restorations are desired, two or more shade tabs can be placed on top of each other to estimate the resulting layered shades.

PACKAGE

1 piece

NOTE

The shade guide can be cleaned with water and mild soap. It should not be sterilized. The use of various chemical disinfectants / sterilants may cause damage to the shade guide.

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3737 West 127th Street, Alsip, IL 60803 U.S.A. TEL: +1-708-597-0900
GC ASIA DENTAL PTE. LTD.