G-CEM LinkForce™ TECHNIQUE GUIDE

CEMENTATION TECHNIQUE FOR VENEERS

TRIAL FIT

Remove the temporary restoration and clean thoroughly.



Check the fit and occlusion. As necessary, use G-CEM Try-In Paste.



Remove the restoration and rinse the paste with water.



▶ PREPARATION OF THE RESTORATION

Prepare the restoration according to the manufacturer's instructions.

▶ Glass ceramics

Etch with hydrofluoric acid.



Rinse and dry.



▶ Hybrid ceramics and Composite

Sandblast*
*For CERASMART,
alternatively apply
hydrofluoric acid.



Blow clean with air syringe. Rinse and dry.









synnge.

▶ PREPARATION OF TOOTH SURFACE



Rinse and dry the prepared tooth.



Select either Selective Etch or Total Etch technique. Rinse and dry.







▶ CEMENTATION



Apply cement directly to the bonding surface of the veneer and/or tooth surface. 10



Immediately seat onto prepared tooth. Maintain moderate pressure.



Remove excess. Excess can be tack cured for 1-2 seconds.



700mW/cm²).

13



CEMENTATION TECHNIQUE FOR POSTS AND CORES

▶ PREPARATION OF THE RESTORATION

Prepare the restoration according to the manufacturer's instructions.










▶ PREPARATION OF POST SPACE












▶ CEMENTATION





Remove excess

cement.





CEMENTATION TECHNIQUE FOR INLAYS, ONLAYS, CROWNS AND BRIDGES

TRIAL FIT

Remove the temporary restoration and clean thoroughly.



Check the fit and occlusion. As necessary, use G-CEM Try-In Paste.



Remove the rinse the paste with



▶ PREPARATION OF THE RESTORATION

Prepare the restoration according to the manufacturer's instructions.

▶ Glass ceramics

Etch with hydrofluoric acid.



Rinse and dry.



▶ Metal, Zirconia, Alumina, Hybrid ceramics and Composite

Sandblast* *For CERASMART, alternatively apply hydrofluoric acid.



Blow clean with air syringe. Rinse and dry.







syringe.





▶ PREPARATION OF TOOTH SURFACE











► Alternative dual-cure mode









▶ CEMENTATION



LinkForce directly into the restoration.



onto prepared tooth/abutment. Maintain moderate pressure.



Remove excess while maintaining moderate pressure. Excess can be tack cured for 1-2 seconds for an easier excess remarks an easier excess remarks and excess remarks remarks and excess remarks remarks rem excess removal





polish margins if

necessary.

