Initial IQ

JUNE 2012

ONE BODY,
LAYERING-OVER-METAL
LAYERING-OVER-ZIRCON

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GC Initial IQ • One Body Layering-over-Metal, Farbtabellen

Colour Charts
Farbtabellen

Choose from the following components to create the most popular shades.
Please find detailed combination recommendations on page 10.

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<i>initial</i> ^{IQ} One Body Layering-over-Metal, Body Powders					
Body Powders	4	A	B	C	D

initial^{IQ} Lustre Pastes / Lasurpasten

GC Initial IQ • One Body Layering-over-Metal, Opaque Layering GC Initial IQ • One Body Layering-over-Metal, Opakerschichtung

Cast or CAD/CAM made metal frameworks are finished using tungsten carbide metal burs or ceramic-bonded grinding instruments and oxidized according to alloy manufacturers instructions.

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Apply the first Opaque layer (wash) with the normal GC Initial MC Paste Opaques in a thin layer using a flat brush. Before firing, fluorescent crystals can be applied on the Opaque surface.

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After firing, these Fluo Crystals are the base for a very important equal layer of the 2nd Opaque firing.

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Fire second Opaque layer as usual, at critical areas the coloured Opaque Modifiers should be used. The Opaque layer should have a slightly shiny surface.

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Application of Body Powder over the entire surface.

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The restoration is fired following the firing parameters of the 1st Body firing (see firing chart). The fired ceramic surface shows a slight gloss.

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GC Initial IQ • Lustre Pastes Application

GC Initial IQ • Lasurpastenanwendung

Before applying the Lustre Coating, the ceramic restoration is shaped and contoured using adjusted diamond burs. After prepolishing with a ceramic polisher the finished crowns and bridges are sand-blasted with 50 microns Alu-Oxid at 1.5 bar pressure.

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The Lustre Paste Neutral is coated on the entire surface of the sandblasted restoration. This coating is clearly thicker than what is known from a normal Glaze firing.

Note: these Lustre Pastes must not come in contact with water.

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The requested Vita Classic Shade is created by a more or less intensive application of the Lustre Pastes "Body A - D". This process, which only takes two minutes, is everything that is necessary to create Vita Classic Shades.

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GC Initial IQ • Lustre Pastes Application

GC Initial IQ • Lasurpastenanwendung

Tip: The desired surface texture or smoothness of the Lustre Paste layer can be obtained by softly vibrating or condensing the restoration.

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Result after 1 firing.

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For a more detailed description of the possibilities of Lustre Pastes NF, please see the appropriate manual.

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

Important Information:

All GC Initial, IQ One Body Layering-over-Metal materials and GC Initial MC materials are fully compatible and can easily be combined for extended options.

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GC Initial IQ GC Initial IQ


GC Initial IQ One Body Layering-over-Metal capabilities.

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GC Initial IQ • One Body Layering-over-Metal, Firing Chart

GC Initial IQ • One Body Layering-over-Metal, Brenntabelle

	Preheating temp. Start-Temperatur	Drying time Trocknungs-Zeit	Raise of temp. Temperatur-Anstieg	Vacuum Vakuum	Final temp. End-Temperatur	Holding time Haltezeit	Appearance Erscheinungsbild
Oxidation Oxidbrand	Refer to the alloy manufacturer's instructions Je nach Angaben des Legierungsherstellers						
1st Paste Opaque Firing* Washbrand* (Pastenopaker)	550°C	6 min	80°C/min	Yes/Ja	940°C	1 min	Shiny Glänzend
2nd Paste Opaque Firing Opakerbrand (Pastenopaker)	550°C	6 min	80°C/min	Yes/Ja	930°C	1 min	Slightly shiny Leicht glänzend
1st Body firing 1. Bodybrand	600°C	4 min	80°C/min	Yes/Ja	900°C	1 min	Slightly shiny Leicht glänzend
2nd Body firing (correction firing) Korrekturbrand	600°C	3 min	80°C/min	Yes/Ja	890°C	1 min	Slightly shiny Leicht glänzend
Glaze Firing Lustre Pastes Glanzbrand Lustre Pastes	480°C	2 min	45°C/min	No/Nein	810°C	1 min	Shiny Glänzend

Note: Ceramic furnaces of different types may vary significantly in real firing temperature compared to indicated temperature. Please periodically calibrate your furnace and do test firings if necessary.


**Please increase temperature of first Opaque firing by +20°C when using Non-Precious Metals.*

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Related Products - INmetalbond

Verwandte Produkte - INmetalbond


Thin layer bonder between alloy and ceramic. The INmetalbond avoids degassing of metal oxides during various firings and neutralizes differences in CTE. The metal framework should be prepared according to the manufacturers instructions.

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Stir before using. Apply the bonder in a thin layer, but completely mask the framework. Important: Do not moisten dried Bonder. Close the cap after using.

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The Bonder has a yellowish, slight shiny surface. Note: The colour of the fired bonder depends on the alloys' compounds. When using non precious alloys with the INmetalbond, the Opaque washbake does not need to be increased by 20°C.

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The above mentioned firing parameters are only guidelines and therefore always need to be adjusted to the firing furnace and its correct functionality. Most important is to obtain the right firing result. These firing parameters can only be used as guidelines.

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GC Initial IQ • One Body Layering-over-Zircon, Farbtabellen

Colour Charts Farbtabellen

Choose from the following components to create the most popular shades.
Please find detailed combination recommendations on page 23.

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<i>initial</i> _{IQ} One Body Layering-over-Zircon, Body Powders						
Body Powders	5	A	B	C	D	Bleach

*initial*_{IQ} Lustre Pastes / Lasurpasten

GC Initial • Zircon - Light Reflective (LR) Liner GC Initial • Zircon - Light Reflective (LR) Liner


Using the Zr-LR Liners, you can easily adjust the shade of different coloured ZrO₂ materials for a systematic approach. For example, to achieve the Classic Vita shades, see the colour combination charts page 23.

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A single-colour application (a) of the Zr-LR Liners is possible, as well as a multi-colour application (b), or an additional coloration with GC Initial IQ Lustre Pastes NF (c).

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
GC Initial • Zircon - Light Reflective (LR) Liner

GC Initial • Zircon - Light Reflective (LR) Liner

Apply Zr-LR Liner evenly thin, but covering.
Refer to the firing chart on page 24.

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Tip: Use more intense Zr-LR Liner to modify the basic shades partially. To achieve maximum modification, we also recommend the use of GC Initial IQ Lustre Pastes on the Liner.

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Apply one Body Powder, building up the entire crowns to anatomical contour (A, B, C, D or BL Powders).

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Tip: *Good condensing before firing reduces shrinkage.*

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GC Initial IQ • One Body Layering-over-Zircon, Full Contour Build Up Procedure

GC Initial IQ • One Body Layering-over-Zircon, Vollanatomische Schichtung

For corrections, use the same powder again. Fire according to firing chart on page 24.

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Correct the shape by appropriately grinding as usual.

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*Sandblast before coloration by Lustre Pastes with Glaze firing
- Alox 50 μ m, 2bar pressure.*

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Finally, apply Lustre Pastes as described on page 8.

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GC Initial IQ • One Body Layering-over-Zircon, Recommendations

GC Initial IQ • One Body Layering-over-Zircon, Hinweise



Highly aesthetic result after Glaze Firing with Lustre Pastes.

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Tip: For optimal aesthetic results, the ZrO_2 framework of the pontic should be extended to basal as far as possible (within the common guidelines of hygienic requirements). If the framework is not appropriately extended, the cervical area of the pontic will have a grayish appearance. Proper framework design, anticipating the transparency, and translucency effects of the veneering ceramics is key for the best aesthetic design.

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

Important Information:


All GC Initial, IQ One Body Layering-over-Zircon materials and GC Initial Zr-FS materials are fully compatible and can easily be combined for extended options.

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GC Initial IQ • One Body Layering-over-Zircon, Firing Chart

GC Initial IQ • One Body Layering-over-Zircon, Brenntabelle

	Preheating temp. Starttemperatur	Drying time Trocknungszeit	Raise of temp. Temperaturanstieg	Vacuum Vakuum	Final temp. Endtemperatur	Holding time Haltezeit
Light Reflective Liner Firing Light Reflective Liner-Brand	450°C	6 min	45°C/min	Yes / Ja	860°C	1 min
1st Body firing 1. Bodybrand	450°C	5 min	50°C/min	Yes / Ja	810°C	1 min
2nd Body firing (correction firing) Korrekturbrand	450°C	4 min	50°C/min	Yes / Ja	805°C	1 min
Glaze Firing Lustre Pastes Glanzbrand Lustre Pastes	480°C	2 min	45°C/min	No / Nein	780°C	1 min

Note: Ceramic furnaces of different types may vary significantly in real firing temperature compared to indicated temperature. Please periodically calibrate your furnace and do test firings if necessary.

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Physical Properties and Shelf Life

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A short overview of the capacity of the GC Initial IQ One Body Layering-over-Zircon System

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Vita® is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.
GC NORDIC AB
GC AUSTRIA GmbH