Master all your challenges

*initial*TM

LF

- EN Technical Manual
- DE Technisches Handbuch
- FR Manuel Technique
- IT Manuale tecnico

AUGUST 2021

GC

EN

DE

Table of Contents	Inhalts-verzeichnis	Page
Introduction & Intended Use	Einleitung & bestimmungsgemaser Gebrauch	4
True-to-nature ...	Naturidentisch ...	5
Colour Combination Chart	Farbtabelle	5
GC Initial LF System Colour Chart	GC Initial LF System Farbkarte	6
Standardized build up procedure Anterior tooth Molar	Standardisierte Schichttechnik Frontzahn Molar	7-11 12-17
GC Initial – Colour combination charts	GC Initial – Farbuordnungstabellen	18-19
Polychromatic build up procedure	Polychrome Schichttechnik	20-28
Firing instructions	Brennanleitung	29
INmetalbond	INmetalbond	30-31
Physical properties & Shelf Life	Physikalische Eigenschaften & Haltbarkeitsdauer	32
Alloy list for GC Initial LF	Legierungsliste für GC Initial LF	33

EN Introduction & Intended Use

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
EN Thank you for choosing our GC Initial LF veneering ceramic.

GC Initial LF is a low-fusing veneering ceramic with EN ISO 6872 for the ceramic veneering of precious and non-precious metal substructures with a coefficient of thermal expansion of the metal substructure in the range of $13,8 - 14,9 \cdot 10^{-6} \text{ K}^{-1}$ (25-500°C). It is also ideal for correcting restorations made with GC Initial MC or as a build up material on GC Initial PC copings. This instruction for use will give you an impression of how easy it is to get a convincing aesthetic result with a minimum workload that highlights the excellent features of this veneering ceramic. The mentioned firing parameters are only guidelines and therefore always need to be adjusted to the firing furnace and its correct functionality. Most important is to obtain the right firing result. These firing parameters can only be used as guidelines.

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Description	No. of Bottles	Colour Code
Cervical Translucent	5	CT-21: Light
		CT-22: Yellow Soft
		CT-23: Pale Orange
		CT-24: Yellow
		CT-25: Brown
Shoulder Transpa	7	ST-30: Neutral White
		ST-31: Cream
		ST-32: Light Yellow
		ST-33: Yellow
		ST-34: Pink/Beige
		ST-35: Olive Brown
		ST-36: Light Brown
Shoulder Opaque	3	SO-37: Ivory
		SO-38: Orange
		SO-39: Brown
INside	11	IN-41: Flamingo
		IN-42: Terracota
		IN-43: Sun
		IN-44: Sand
		IN-45: Havanna
		IN-46: Brasil
		IN-47: Sienna
		IN-48: Kurkuma
		IN-49: Maracuja
		IN-50: Curry
		IN-51: Olive
Glaze	1	GL
Correction Powder	1	COR
Bleach	4	BLD-1: Light
		BLD-2: White
		BLD-3: Xwhite
		BL-E: Bleach enamel
Gum	1	GU: Gum Universal

EN Standardized build up procedure
DE Standardisierte Schichttechnik

EN Cast metal frameworks are finished using tungsten carbide metal burs or ceramic-bonded grinding instruments and oxidized according to alloy manufacturer's instructions. INmetalbond application, see page 30.

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
EN Fired 2nd Opaque layer with Opaque Modifiers, the Opaque layer should have a slightly shiny surface.

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EN Standardized build up procedure
DE Standardisierte Schichttechnik


EN Inside, Primary Dentin, is applied on the cervical part and Opaqus Dentin on the incisal border of the framework. (see reference chart)

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EN Dentin material in its correct shade is applied.

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EN The complete Dentin layering in the desired inner structure.

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EN Enamel (E) and / or Translucent (T) are applied on top of the Clear Fluorescence (CL-F) up to the final shape.

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

EN Standardized build up procedure DE Standardisierte Schichttechnik


EN The marginal ridges are lined with the same Dentin mixture.

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EN The palatal fossa of the restoration is covered with the high chroma INside material.

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EN Final oral build up with the same Enamel (E) mixture that has been used on the labial side.

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EN Standardized build up procedure
 DE Standardisierte Schichttechnik

EN The restoration is fixed following the firing parameters of the 1st Dentin firing (see firing instructions). The fired ceramic surface shows a slight gloss.

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EN The proximal areas are supplemented with the same material used for the Dentin / Enamel layering. The restoration is finally adjusted using Enamel Transparent materials. Cervical Translucent (CT) can be used as an alternative for the Dentin / Enamel layering.

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EN Completed "GC Initial" restoration after glazing. Application of individualized characterizations with ceramic stains* is possible.

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EN Remove the oxidation layer with disposable, pure 110 μ Al₂O₃ only and with a maximum pressure of 2 bar. Always observe the instructions for use of the corresponding alloy manufacturer. INmetalbond application, see page 30.

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
EN Standardized build up procedure
DE Standardisierte Schichttechnik


EN The margins and the occlusal fissures line are covered with the selected INside colour (see system colour chart).

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
EN Opaqus Dentin (OD) or Fluo Dentin, FD-91 – FD-93, is applied on the occlusal light refraction areas (see system colour chart).

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EN Standardized build up procedure
DE Standardisierte Schichttechnik

EN Standard Dentin powder is applied in its final form.

The Dentin layering is built up in relatively thin layers due to the fact that these materials all possess a higher chroma level.

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
- EN Standardized build up procedure
- DE Standardisierte Schichttechnik


EN This Dentin base is now covered with Clear Fluorescence material, CL-F. This layer is 0,2 mm in depth and is a copy of human dentition.

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

EN Enamel (E) and / or Translucent (T) are applied on top of the Clear Fluorescence layer up to the occlusal contacts of the antagonists.

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EN Standardized build up procedure
DE Standardisierte Schichttechnik

EN The occlusal plane of the articulator should be raised slightly, (approx. 0,2 mm) in order to build up the final morphology with EO-15.

The use of this versatile special Enamel is crucial to creating depth in the build up, because of it's ability to imitate the natural reflection of the Dentin in the incisal and occlusal edges.

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EN GC Initial INside – Field of application and reference chart

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EN High Fluorescence Dentin

Field of application and reference chart

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EN The Shoulder Opaque SO-37 should be used as a first basic layering at the border frame / tooth shoulder. The percentages used in this list can vary depending on the thickness of the total ceramic shoulder. **BASIC RULE:** The thinner the total ceramic shoulder, the more high chroma shoulder powders should be used. Furthermore, all shoulder powders can be mixed with the INside powders to intensify or to individualize the shoulder shades until a maximum of 20% INside powders of this total mixture. In this case the final temperature of shoulder firing should be reduced by 10°C to 20°C, depending on the amount of INside powders mixed with the shoulder porcelain.

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EN Polychromatic build up procedure
 DE Polychrome Schichttechnik

EN Reduced metal frame on the shoulder for optimal light transmission at the margins.
 INmetalbond application, see page 30.

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EN Polychromatic build up procedure
DE Polychrome Schichttechnik

EN First ceramic shoulder build up with a thin layer of high fluorescence Opaque shoulder material.

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EN 2nd shoulder firing, towards the margin line, is carried out with a more translucent shoulder porcelain, ST, to allow light transmission through the gingiva.

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
EN Polychromatic build up procedure

DE Polychrome Schichttechnik


EN The ceramic shoulder after firing.

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EN Opaqus Dentin (OD) or Fluo Dentin (FD-91 – FD-93) is applied on the incisal border of the coping towards the incisal edge.

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EN The margins and the body part of the crown are covered with the selected INside colours (see colour combination chart)

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EN Polychromatic build up procedure
DE Polychrome Schichttechnik


EN Standard Dentin powder is mixed with 10% of INside powder (see colour combination chart) and very thinly applied (about 50% of the usual thickness).

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
EN Dentin is applied in its final form. This way of working allows a true-to-nature build up procedure of the Transpa / Enamel layering.

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
EN Polychromatic build up procedure

DE Polychrome Schichttechnik


EN By means of this Transpa / Enamel layering, the incisal, approximal end form has been created.

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EN By means of a "cut-back", the incisal part is reduced up to the Opacus Dentin layering or Fluo Dentin layering.

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EN The reduced area is reconstructed with FD-91 in the desired Mamelon structure. This is a very important point for colour synchronization of different base structures placed next to each another.

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
EN Polychromatic build up procedure
DE Polychrome Schichttechnik

EN Displaying the colour intensity of the Mamelon structure is realized with ceramic stains* (see GC Initial LF system colour chart).

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
* GC Initial Spectrum Stains / GC Initial INvivo/INSitu Stains

EN This base is now covered with clear, fluorescence material, CL-F. This layer is 0,2 mm in depth and is a copy of human dentition.

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
EN The Enamel layering starts with applying blue or grey Enamel Opal (EOP3 / EOP4) in the incisal and / or approximal area.

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
EN Polychromatic build up procedure

DE Polychrome Schichttechnik


EN The Enamel / Transpa layering, where special incisal and transparent are alternatively placed next to another (e.g. E58, TO, EOP2, CL-F)

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EN The final contouring is made out of a mixture of Enamel Opalescence (EOP) and standard Enamel (E). Depending on the age of the patient, the translucent level of the Enamel layer can be modified with TN, TO & CL-F. To create the milky white incisal contour, EO-15 should be used.

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EN Polychromatic build up procedure
DE Polychrome Schichttechnik

EN The palatal fossa of the restoration is covered with a high-chroma INside colour (see colour combination chart).

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EN Polychromatic build up procedure

DE Polychrome Schichttechnik

EN The final functional contouring of the oral side is done with Enamel Occlusal (EO-15).

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EN The result of an always returning layering system.

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EN * When firing on NPA (non-precious alloys) increase final temperature with 20°C.

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
EN The above mentioned firing parameters are only guidelines and therefore always need to be adjusted to the firing furnace and its correct functionality. Most important is to obtain the right firing result. These firing parameters can only be used as guidelines.

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EN Thin layer bonder between alloy and ceramic. The INmetalbond avoids degassing of metal oxides during various firings and neutralizes differences in CTE. The metal framework should be prepared according to the manufacturer's instructions.

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EN Stir before using. Apply the bonder in a thin layer but mask the framework complete. Important: Do not moisten dried Bonder. Close the cap after using.

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EN Related Products – INmetalbond
 DE Vergleichbare Produkte – INmetalbond

EN The Bonder has a yellowish, slight shiny surface.
 Note: The colour of the fired bonder depends on the alloys compounds. When using non precious alloys with the INmetalbond, the Opaque washbake does not need to be increased by 20°C.

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
EN GC Initial LF is a low-fusing veneering ceramic for all types of precious and non-precious alloys. Make sure that the CTE of the alloy you use for casting or milling the substructure stays strictly within the range of 13,8 - 14,9 · 10⁻⁶ K⁻¹ (25-500 °C). Calibrate your furnace according to the manufacturer's instruction in order to obtain the best results and to be able to make optimal use of the properties. As each dental ceramic is sensitive to contamination, take care that you work in a clean environment.

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EN Medical Device Regulation Requirement

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

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Notes

Notes

Lined writing area with multiple horizontal lines.

 EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240
Interleuvenlaan 33, B-3001 Leuven
Tel. +32.16.74.10.00, Fax. +32.16.40.02.14
info.gce@gc.dental, <https://europe.gc.dental>

RESPONSIBLE MANUFACTURER IN CANADA
GC AMERICA INC.
3737 W. 127th Street, Alsip, IL 60803 U.S.A.

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
Tel: +1-708-597-0900
www.gcamerica.com

GC Germany GmbH
Seifgrundstraße 2, D-61348 Bad Homburg
Tel. +49.61.72.99.59.60, Fax. +49.61.72.99.59.66.6
info.germany@gc.dental, <https://europe.gc.dental/de-DE>

GC ITALIA S.r.l.
Via Calabria 1, I-20098 San Giuliano Milanese
Tel. +39.02.98.28.20.68, Fax. +39.02.98.28.21.00
info.italy@gc.dental, <https://europe.gc.dental/it-IT>

GC UNITED KINGDOM Ltd.
Coopers Court, Coopers Court, Newport Pagnell, UK-Bucks. MK16 8JS
Tel. +44.1908.218.999, Fax. +44.1908.218.900
info.uk@gc.dental, <https://europe.gc.dental/en-GB>

GC FRANCE s.a.s.
8 rue Benjamin Franklin, F-94370 Sucy en Brie Cedex
Tel. +33.1.49.80.37.91, Fax. +33.1.45.76.32.68
info.france@gc.dental, <https://europe.gc.dental/fr-FR>

GC IBÉRICA - Dental Products, S.L.
Edificio Codesa 2
Playa de las Américas, 2, 1º, Of. 4, ES-28290 Las Rozas, Madrid
Tel. +34.916.364.340, Fax. +34.916.364.341
comercial.spain@gc.dental, <https://europe.gc.dental/es-ES>

GC AUSTRIA GmbH
Tallak 124, A-8103 Gratwein-Strassengel
Tel. +43.3124.54020, Fax. +43.3124.54020.40
info.austria@gc.dental, <https://europe.gc.dental/de-AT>

GC Europe NV - Benelux Sales Department
Researchpark Haasrode-Leuven 1240
Interleuvenlaan 33, B-3001 Leuven
Tel. +32.16 74.18.60
info.benelux@gc.dental, <https://europe.gc.dental/fr-BE>

GC EUROPE N.V. - East European Office
Siget 19b, HR-10020 Zagreb
Tel. +385.1.46.78.474, Fax. +385.1.46.78.473
info.eeo@gc.dental, www.eeo.gceurope.com

GC NORDIC AB - Finnish Branch
Lemminkäisenkatu 46, FIN-20520 Turku
Tel. +358.40.900.07.57
info.finland@gc.dental, <https://europe.gc.dental/fi-FI>

GC NORDIC AB - GC Nordic Danish Branch
Scandinavian Trade Building
Gydevang 34-41, DK-3450 Allerød
Tel. +45 51 15 03 82
info.denmark@gc.dental, <https://europe.gc.dental/da-DK>

GC NORDIC AB
c/o Lundin Revisionbyrå
Erik Dahlbergsgatan 11B, SE-411 26 Göteborg
Tel. +46.768.54.43.50
info.nordic@gc.dental, <https://europe.gc.dental/sv-SE>

GC AUSTRIA GmbH - Swiss Office
Zürichstrasse 31, CH-6004 Luzern
Tel. +41.41.520.01.78, Fax. +41.41.520.01.77
info.switzerland@gc.dental, <https://europe.gc.dental/de-CH>

GC AUSTRALASIA DENTAL PTY LTD
1753 Botany Rd, Banksmeadow NSW 2019, Australia
Tel: +61-2-9301 8200, Fax: +61-2-9316 4196

GC SOUTH AMERICA
Rua Heliadora, 399, Santana - São Paulo, SP, BRASIL
CEP: 02022-051 - TEL: +55-11-2925-0965 - CNPJ: 08.279.999/0001-61
RESP. TÊC: Erick de Lima - CRO/SP 100.866

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529541
Tel: +65 6546 7588

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EN V-Shade is referring to the Vita® Shade Guide. Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.