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1. 大島克郎, 竹井利香, 安藤雄一. 就業歯科技工士数の将来推計. 日歯医療管理会誌 2019;54(3):199-207.
2. 厚生労働省: 歯科用貴金属随時改定について. 中医学 総-6資料 2024. <https://www.mhlw.go.jp/content/12404000/001209362.pdf>.
3. Govare N, Contrepolis M. Endocrowns: A systematic review. J Prosthet Dent 2020;123(3):411-418.
4. 公益社団法人日本補綴歯科学会: 保険診療における CAD/CAM冠の診療指針 2024.
5. Belleflamme MM, Geerts SO, Louwette MM, Grenade CF, Vanheusden AJ, Mainjot AK. No post-no core approach to restore severely damaged posterior teeth: An up to 10-year retrospective study of documented endocrown cases. J Dent. 2017 ;63:1-7.
6. Fages M, Bennisar B. The endocrown: a different type of all-ceramic reconstruction for molars. J Can Dent Assoc. 2013;79:d140.
7. Demachkia AM, Velho HC, Valandro LF, Dimashkieh MR, Samran A, Tribst JPM, de Melo RM. Endocrown restorations in premolars: influence of remaining axial walls of tooth structure and restorative materials on fatigue resistance. Clin Oral Investig. 2023;27(6):2957-2968.
8. Jalali S, Asgari N, Pirooz P, Younespour S, Atri F. Comparison of Clinical Efficacy of CAD/CAM Endocrowns Made of Feldspathic, Zirconia Lithium Silicate, and Lithium Disilicate: A Two-year Mixed Cohort Study. J Dent. 2024;24:105019.
9. 松村秀雄, 小峰太. 支台歯形成. 矢谷博文, 三浦宏之, 細川隆司, ほか編, クラウンブリッジ補綴学 第6版, 東京: 医歯薬出版; 2021, 114-125.
10. Zheng Z, Sun J, Jiang L, Wu Y, He J, Ruan W, Yan W. Influence of margin design and restorative material on the stress distribution of endocrowns: a 3D finite element analysis. BMC Oral Health. 2022;22(1):30.
11. Pedrollo Lise D, Van Ende A, De Munck J, Umeda Suzuki TY, Cardoso Vieira LC, Van Meerbeek B. Biomechanical behavior of endodontically treated premolars using different preparation designs and CAD/CAM materials. J Dent. 2017;59:54-61.
12. Dartora NR, de Conto Ferreira MB, Moris ICM, Brazão EH, Spazin AO, Sousa-Neto MD, Silva-Sousa YT, Gomes EA. Effect of Intracoronary Depth of Teeth Restored with Endocrowns on Fracture Resistance: In Vitro and 3-dimensional Finite Element Analysis. J Endod. 2018;44(7):1179-1185.
13. Soliman M, Alzahrani G, Alabdualataif F, Eldwakhly E, Alsamady S, Aldegheishem A, Abdelhafeez MM. Impact of Ceramic Material and Preparation Design on Marginal Fit of Endocrown Restorations. Materials (Basel). 2022;15(16):5592.
14. Shin Y, Park S, Park JW, Kim KM, Park YB, Roh BD. Evaluation of the marginal and internal discrepancies of CAD-CAM endocrowns with different cavity depths: An in vitro study. J Prosthet Dent. 2017;117(1):109-115.
15. Hiraba H, Nishio K, Takeuchi Y, Ito T, Yamamori T, Kamimoto A. Application of one-piece endodontic crowns fabricated with CAD-CAM system to molars. Jpn Dent Sci Rev. 2024;60:81-94.
16. Tribst JPM, Lo Giudice R, Dos Santos AFC, Borges ALS, Silva-Concilio LR, Amaral M, Lo Giudice G. Lithium Disilicate Ceramic Endocrown Biomechanical Response According to Different Pulp Chamber Extension Angles and Filling Materials. Materials (Basel). 2021;14(5):1307.
17. Zhu J, Wang D, Rong Q, Qian J, Wang X. Effect of central retainer shape and abduction angle during preparation of teeth on dentin and cement layer stress distributions in endocrown-restored mandibular molars. Dent Mater J. 2020;39(3):464-470.
18. Yamauchi S, Shipper G, Buttke T, Yamauchi M, Trope M. Effect of orifice plugs on periapical inflammation in dogs. J Endod. 2006;32(6):524-6.
19. Ikemoto S, Komagata Y, Yoshii S, Masaki C, Hosokawa R, Ikeda H. Impact of CAD/CAM Material Thickness and Translucency on the Polymerization of Dual-Cure Resin Cement in Endocrowns. Polymers (Basel) 2024;16(5).
20. Thomas RM, Kelly A, Tagiyeva N, Kanagasingam S. Comparing endocrown restorations on permanent molars and premolars: a systematic review and meta-analysis. Br Dent J. 2020.
21. Al-Dabbagh RA. Survival and success of endocrowns: A systematic review and meta-analysis. J Prosthet Dent 2021;125(3):415 e1- e9.
22. European Society of Endodontology developed by: Mannocci F, Bhuva B, Roig M, Zarow M, Bitter K. European Society of Endodontology position statement: The restoration of root filled teeth. Int Endod J. 2021 Nov;54(11):1974-1981.
23. 駒形裕也, 正木千尋, 吉居慎二, 池田弘, 渡辺崇文, 畑健太郎, 赤間廣輔, 加来伸哉. 低侵襲な修復治療を実現するエンドクラウン. デンタルダイヤモンド 2023; 48 (12), 27-47.

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