

initial™ Zirconia Coloring Liquid

Customised, your way!

Faster, more predictable
and highly individualised
following the Initial IQ philosophy

Smaller inventory, greater gains
Obtain every shade fast and easy
with a reduced inventory



Brush or dip
Two ways, your
call



Since 1921
100 years of Quality In Dental



Set the tone for efficiency - following the GC Initial IQ philosophy.

The **Initial Zirconia Coloring Liquid** system consists of 12 water-based colouring liquids following the 'IQ philosophy'.

It is optimised for use with the **Initial Zirconia Disk HT**, a universal white zirconia that offers an ideal balance between translucency, strength and aesthetic possibilities.

The Initial Zirconia Coloring Liquid system guarantees the replication of natural tooth colours and effects prior to sintering while simplifying inventory management. Simply apply the liquids onto the monolithic or anatomically reduced restoration after milling using the painting or dipping technique.

This framework in the right colour tone can then be finished using the **Initial IQ ONE SQIN** concept, ensuring a fast, aesthetically superior ceramic finish through painting and micro-layering.



Paint it or dip it - Two ways to accomplish

Ultimate frame customisation

Precisely controlled and localised colour penetration
Reduces the need for post-sintering colour adjustments.




Deep, even colour distribution

Consistent and reproducible results in a fast and easy manner.



Reduced amount of liquids

- 4 cervical shades
- 4 body shades
- 4 effect shades



Cervical A


Cervical B

Cervical C

Cervical D

Easy diluting chart

to achieve all 16 V-shades.



Body A

Body B

Body C

Body D

Provide a **subtle and diffuse undertone** - a perfect base for the IQ ONE SQIN technique.



Effect Blue

Effect Grey

Effect Neutral

Effect Pink

4 marking colours

facilitate the visualisation during application.





Benefits at a glance

- For complete or partial characterisation of zirconia frameworks, either by **painting** or **dipping**
- **Nature analogue** colour tones to limit the required finishing time
- **Enhanced aesthetics** by customisation of the pre-sintered zirconia structure
- **Faster**, more predictable and highly individualised results

Perfect synergies



Initial Zirconia Disk



Initial Zirconia Coloring Liquids



Initial IQ ONE SQIN