

# G-ænial™ Universal Injectable

Universal Light-Cured Radiopaque Injectable Composite

Transform the way  
you work




Since 1921

100 years of Quality in Dental



An injectable high-strength ultra-fine particle composite with ideal viscosity, handling and adaptation characteristics that may be used for long-lasting esthetic restorations.

The homogeneously dispersed ultra-fine barium-glass particles formula (150nm particle size, 69% wt filler loading) provides high flexural strength and wear resistance for a durable and long lasting restoration.



Metal tip may be bent with fingers if necessary to better access posterior distal region.

Cavity Depth	Product Recommendation
Up to 2 mm	<b>G-aenial™ Flo X</b> or
	<b>G-aenial™ Universal Injectable</b>
Up to 4-6+ mm	<b>G-aenial™ BULK Injectable with G-aenial Sculpt®</b> (veneer layer optional)

Start injecting with our strongest direct restorative ever\*. G-aenial™ Universal Injectable features high flexural strength, for protection against chipping and fracture, and outstanding wear resistance for excellent durability and gloss retention.



#### G-aenial™ Universal Injectable Syringe Refills

Contains: Two 1.7g (1.0mL) syringes, 10 long dispensing tips and 10 short dispensing tips.

- |            |            |            |             |
|------------|------------|------------|-------------|
| 012364 A1  | 012365 A2  | 012366 A3  | 012367 A3.5 |
| 012368 A4  | 012369 B1  | 012370 B2  | 012371 CV   |
| 012372 CVD | 012373 AO1 | 012374 AO2 | 012375 AO3  |
| 012376 JE  | 012377 AE  | 012378 XBW | 012379 BW   |

#### G-aenial™ Universal Injectable Unitip Refills

Contains: 15 0.16mL unitips.

- |            |            |            |             |
|------------|------------|------------|-------------|
| 012348 A1  | 012349 A2  | 012350 A3  | 012351 A3.5 |
| 012352 A4  | 012353 B1  | 012354 B2  | 012355 CV   |
| 012356 CVD | 012357 AO1 | 012358 AO2 | 012359 AO3  |
| 012360 JE  | 012361 AE  | 012362 XBW | 012363 BW   |

425000 Unitip APPLIER II

Scan for a Free Lunch & Learn!

