



Catalog N°6

Facebows, Lab, Miscellaneous



Experts in Orthodontics
GC Orthodontics Europe GmbH
www.gcorthodontics.eu

Facebows, Lab, Miscellaneous

- DentaSonic.....219
- Interproximal reduction.....219
- Laboratory222
 - Smile Rite Thermoforming.....222
 - Extraoral Traction.....222
- Miscellaneous223
 - Band Pushers223
 - Cheek Retractor.....223
 - Cheek Retractor with vaccuming system NOLA.....224
 - Protective Wax224
 - Protective Silicon224
 - Mirrors225
 - Lip Protection.....225
 - Elastics Placement.....225
 - Elasti-bites.....225



Interproximal reduction

For safe, precise surface reduction with minimal enamel damage, our professional materials for MI-IPR, reduce inter-proximal spaces with minimal deterioration.



Profi Kit 2

Profi Kit 2 Box includes single and double-sided DiaStrips, for inter-proximal space reduction, proximal finishing and polishing in orthodontics, aligners and restorative dentistry. Profi Kit 2 also features a new water-cooled oscillating handpiece with a 4:1 ratio and a second head with an alternative movement.

Description	P/N
4 DiaStrips, 2 sides • 4 DiaStrips 1 side • 2 Polishing PoliStrip • 2 polishing burs • 1 water cooled hand piece, 4:1 • 1 tête oscillante • 1 measuring gauge set: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5mm • 1 extracting tip • 1 Spray	55-BD02-0000



Basic set nwc (no water cooled)

Basic reduction set of 4:1.

Description	P/N
Reduction: 4:1 (not cooled) • 1 assortment of 4 stripes • 1 measuring gauge set: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5mm • 1 extracting tip	55-BD04-0000



Double-sided DiaStrip

Double-sided DiaStrip oscillating diamond files, allows for inter-proximal reduction on both sides at the same time. Used for proximal finishes and orthodontic polishing.

Double-Sided	DiaStrip	DiaStrip	DiaStrip	DiaStrip
P/N 5/pk	55-FDM0-5215	55-FDM0-5225	55-FDM0-5240	55-FDM0-5260
P/N 12/pk	55-FDM0-1215	55-FDM0-1225	55-FDM0-1240	55-FDM0-1260
Thickness	0.1mm	0.15mm	0.2mm	0.3mm



Single-sided DiaStrip

Single-sided DiaStrip oscillating diamond files, allows inter-proximal reduction on one side only. Used for proximal finishes and orthodontic or restorative polishing.

Simple-Sided	DiaStrip	DiaStrip	DiaStrip	DiaStrip
P/N 5/pk	55-FDM0-5115	55-FDM0-5125	55-FDM0-5140	55-FDM0-5160
Thickness	0.1mm	0.15mm	0.2mm	0.3mm



DiaSaw

Diamond oscillating DiaSaw files. Opens contact points safely and quickly (without soft tissue damage) and reduces tooth size in orthodontics.

Double sides	DiaSaw	DiaSaw	DiaSaw
P/N 5/pk	55-SDM0-5315	55-SDM0-5325	55-SDM0-5340
Thickness = +/- 0.03mm	0.1mm	0.15mm	0.2mm



Set



Double-Sided	DiaStrip Set 5	DiaSaw Set 5	Aligner Set
P/N	55-FDM0-5200	55-SDM0-5300	55-FDA0-5300
Grain µm:	1x1,0.15,0.2,0.3mm +1PoliStrip	1x1,0.15,0.2mm, +2 PoliStrip	1x1,0.15,0.2mm, +2 PoliStrip
Coating side:	2 sides	2 sides	2 sides
Refill pack:	5 pieces	5 pieces	5 pieces

PoliStrip

Made from flexible resin, these strips are for ultimate proximal polishing on both sides in just 1 passage.

PoliStrips are sold in blister packs of 12 (2x6) with 6 0.10 mm strips and 6 0.15 mm strips.

55-FDRA-1200



DS-Debonding Polisher

Easily and safely remove orthodontic adhesive residue following bracket debonding. Due to the hardness and abrasiveness of materials used with a rotary bur or a metal instrument, there is always a risk when debonding of damaging the natural tooth structure. Our DS-Debonding Polisher features micro particles deposited on a synthetic rubber matrix that safely removes adhesive residue in minutes.

55-D000-6B00



Measuring Gauge

Thickness: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5 mm

55-G000-1000



Inox SteryStrip

Stainless steel sterilization set for 6 polishers.

55-TDM0-6000

Laboratory



Smile Rite Thermoforming


All Smile Rite Thermoforming plastic sheets, regardless of thickness, are pre-dried and packaged in aluminium bags of 25 pieces. The water absorption by the material has been eliminated as well as the bubble problem during thermoforming has been resolved.

Hard	Rond 120 mm	Rond 125 mm	Square 125mmX125mm
1,5 mm (x 50)	—	92-1015-R125	92-1015-S125
1,0 mm (x100)	92-1010-R120	92-1010-R125	92-1010-S125
0,80 mm (x100)	—	92-1008-R125	—

Soft	Rond 120 mm	Rond 125 mm	Square 125 mm X125mm
1,0 mm (X100)	92-1110-R120	92-1110-R125	92-1110-S125

Extraoral Traction

These traction devices are made from washable hypoallergenic materials. Comfortable, in vivid colors, they are easily accepted by young patients.



Neck Strap with Safety Release

Description	P/N
Blue/Pink 450 g	55-4450-0029
Blue/Yellow 450 g	55-4450-0028
Black 600 g	55-4600-0027



Head Cap with Safety Release

Description	P/N
Head Cap Denim 450 g	55-5450-0030
Head Cap Color 600 g	55-5600-0035
Head Cap Denim 600 g	55-5600-0030

Miscellaneous



Band Pushers

Description	P/N
Bite Stick, Autoclavable, Orange	50-8020-0000
Bite Stick, Autoclavable, Blue, Lower Tip	50-8022-0000
Bite Stick, Autoclavable (High Temperature), Blue	50-8023-0000



Cheek Retractor

Description	P/N
Opac, Automatic Cheek Retractor, Full View With Extension	55-0240-0012



Description	P/N
Disposable Cheek Retractor, Medium (x 10)	55-0243-0002



Cheek Retractor with vaccuming system NOLA

Description	P/N
Large, white	55-4100-0000
Narrow, red	55-4101-0000
Tongue pusher large by 2	55-4102-0000
Tongue pusher narrow by 2	55-4103-0000
Rétracteur large, white	55-4104-0000
Pipes with connector "Y" white by 4	55-4105-0000
Connector «Y» by 4	55-4106-0000
Silicon pipes par 10	55-4107-0000



Protective Wax

Description	P/N
Protective Wax, Neutral Flavored (x 100)	55-0111-0033



Protective Silicon

Description	P/N
Protective Silicon (x 50)	55-0121-0000



Mirrors

Description	P/N
Set of 4 Premium Mirrors	55-0170-0000
Premium Occlusal Mirror	55-0171-0000



Lip Protection

Description	P/N
Lip Protector (x 10)	55-0130-0000



Elastics Placement

Description	P/N
Instruments for Elastics, Assorted Colors (x 100)	55-0140-0021



Elasti-bites

Description	P/N
Elasti-bites (x 25)	55-0150-0013

