

Catalog N°6

Bands



Bagues


Band Coding.....	157
Sure-Snap Molar Bands.....	160
A-Fit Molar Bands.....	161
Sure-Snap Bicuspid Bands.....	162

Bands



Two types of high-quality molar bands are available. Both bands are manufactured by TOMY® in Japan.

Band Coding

Sure-Snap Molar Bands: How are product numbers (P/Ns) created?




A-Fit Molar Bands: How are product numbers (P/Ns) created?





Sure-Snap Bicuspid Bands: How are product numbers (P/Ns) created?



Sure-Snap Molar Bands

Sure-Snap Molar Bands are made from a special alloy creating a band that offers at the same time solidity and elasticity. The difference in thickness between occlusal and gingival edges creates bands that are:

- Resistant to the pressure of the band-pusher (occlusal edge)
- Elastic enough to fit over the contour of the tooth without deformation (gingival edge)

Sure-Snap features a unique anatomical form for each type of molar, an easy-to-read permanent identification and an electroeroded interior for increased bond strength.

These easy-to-fit bands are available in 32 sizes for first molars and 30 sizes for second molars.

Sure-Snap First Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2									
Intro	1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	3	3	3	3	2	2	2	2	2	1	1	1	1	1	1	1	1
Standard	1	1	1	1	2	2	2	3	4	4	4	5	5	6	6	6	6	6	5	5	4	3	2	2	2	1	1	1	1	1	1	1
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	14	12	10	8	6	4	4	2	2	2	2	1	1	1

Sure-Snap Second Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2								
Intro	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	3	3	3	3	2	2	2	2	2	2	1	1	1	1	1	1
Standard	1	1	1	2	2	2	2	3	4	4	4	5	5	6	6	6	6	6	5	5	4	3	2	2	2	2	1	1	1	1	1
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	16	12	10	8	6	4	4	2	2	2	2	2	1

Upper First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes
Intro Kit: 100 Assorted Bands, 32 Sizes
Standard Kit: 200 Assorted Bands, 32 Sizes
Master Kit: 400 Assorted Bands, 32 Sizes

Upper Second Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes
Intro Kit: 100 Assorted Bands, 30 Sizes
Standard Kit: 200 Assorted Bands, 30 Sizes
Master Kit: 400 Assorted Bands, 30 Sizes

Lower First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes
Intro Kit: 100 Assorted Bands, 32 Sizes
Standard Kit: 200 Assorted Bands, 32 Sizes
Master Kit: 400 Assorted Bands, 32 Sizes

Lower Second Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes
Intro Kit: 100 Assorted Bands, 30 Sizes
Standard Kit: 200 Assorted Bands, 30 Sizes
Master Kit: 400 Assorted Bands, 30 Sizes



A-Fit Molar Bands

The A-Fit band features all the qualities of Sure-Snap. The difference is a higher crown height and a larger choice of sizes in the higher numbers. See page 109 for A-Fit Bands identification.

A-Fit First Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)

Size	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38		
Starter									1	2	2	2	2	2	2	2	2	2	2	2	2															
Intro	1	1	1	1	1	1	2	2	2	2	2	3	3	3	2	2	2	2	2	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	
Standard	1	1	1	2	4	4	4	5	6	6	6	6	6	6	5	5	4	4	3	3	2	2	2	2	1	1	1	1	1	1	1	1	1	1		
Master	2	2	2	6	8	8	8	10	12	12	12	12	12	12	10	10	8	8	6	6	4	4	4	2	2	2	2	2	2	2	2	2	2	2		

A-Fit Second Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
Starter												1	2	2	2	2	2	2	2	2	2	2	2											
Intro	1	1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	3	3	3	3	2	2	2	2	1	1	1	1	1	1	1	1	
Standard	1	1	1	1	2	2	2	3	4	4	4	5	5	6	6	6	6	6	6	5	5	4	3	2	2	2	1	1	1	1	1	1	1	
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	14	12	10	8	6	4	4	2	2	2	2	1	1	1	1	

Upper First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes
 Intro Kit: 100 Assorted Bands, 34 Sizes
 Standard Kit: 200 Assorted Bands, 34 Sizes
 Master Kit: 400 Assorted Bands, 34 Sizes

Lower First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes
 Intro Kit: 100 Assorted Bands, 34 Sizes
 Standard Kit: 200 Assorted Bands, 34 Sizes
 Master Kit: 400 Assorted Bands, 34 Sizes

Upper Second Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes
 Intro Kit: 100 Assorted Bands, 32 Sizes
 Standard Kit: 200 Assorted Bands, 32 Sizes
 Master Kit: 400 Assorted Bands, 32 Sizes

Lower Second Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes
 Intro Kit: 100 Assorted Bands, 32 Sizes
 Standard Kit: 200 Assorted Bands, 32 Sizes
 Master Kit: 400 Assorted Bands, 32 Sizes

Sure-Snap Bicuspid Bands

Manufacturing quality is identical for these bands and for our molar bands.

Band seating is done through a rocking movement, alternating mesial and distal pressures. Available in three anatomical forms:

- First Maxillary Bicuspids (32 Sizes)
- First Mandibular Bicuspids (25 Sizes)
- Second Mandibular Bicuspids (27 Sizes)

Thanks to their easy seating, you will not hesitate to use Sure-Snap Bicuspid Bands instead of brackets for severe rotations or in case of bond failures.

Sure-Snap Upper Bicuspid Kits – Distribution

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter						2	2	2	2	2	2	2	2	2	2	4	4	4	2	2	2	2	2	2	2	2	2					
Intro	2	2	2	2	2	2	2	2	2	2	4	4	4	4	6	6	6	6	6	4	4	4	4	2	2	2	2	2	2	2	2	2
Standard	2	2	2	2	2	4	6	6	6	6	6	8	8	10	12	14	14	14	12	10	8	8	6	6	6	4	4	4	2	2	2	2
Master	4	4	4	6	6	10	10	14	14	16	16	16	18	18	20	22	22	22	20	18	16	16	14	14	12	12	8	8	8	4	4	4

Sure-Snap First Lower Bicuspids Kits - Distribution

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	
Starter			2	4	4	4	4	4	4	4	4	4	4	4	4											
Intro	2	2	2	2	2	2	4	4	4	4	4	6	6	8	8	8	8	4	4	4	4	2	2	2	2	
Standard	2	2	2	2	4	4	8	8	8	8	10	12	12	18	18	18	18	10	8	8	8	4	4	2	2	
Master	4	4	4	6	6	8	8	14	14	18	18	24	24	36	36	36	36	24	24	16	16	8	8	4	4	

Sure-Snap Second Lower Bicuspids Kits - Distribution

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter						1	1	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	1	1
Intro						2	2	2	2	2	2	2	4	4	4	6	6	8	8	8	6	6	4	4	4	2	2	2	2	2	2	2
Standard						2	2	2	4	4	6	6	6	8	10	12	14	16	16	16	14	12	10	8	6	6	6	4	4	2	2	2
Master						2	4	4	6	8	10	12	12	16	20	26	28	34	34	34	28	26	20	16	14	12	10	8	6	4	4	2

Upper Bicuspids, Bands Only

Starter Kit: 50 Assorted Bands, 22 Sizes

Intro Kit: 100 Assorted Bands, 32 Sizes

Standard Kit: 200 Assorted Bands, 32 Sizes

Master Kit: 400 Assorted Bands, 32 Sizes

First Lower Bicuspids, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes

Intro Kit: 100 Assorted Bands, 25 Sizes

Standard Kit: 200 Assorted Bands, 25 Sizes

Master Kit: 400 Assorted Bands, 25 Sizes

Second Lower Bicuspids, Bands Only

Starter Kit: 50 Assorted Bands, 27 Sizes

Intro Kit: 100 Assorted Bands, 27 Sizes

Standard Kit: 200 Assorted Bands, 27 Sizes

Master Kit: 400 Assorted Bands, 27 Sizes