

Catalog N°6

Elastomerics



Experts in Orthodontics
GC Orthodontics Europe GmbH
www.gcorthodontics.eu























Elastomerics

Elastomeric Color Chains.....	199
Orchestra Intra and.....	200
Extraoral Elastics.....	200
GC Ties Elastomeric Ligatures	201
Hollow or Solid Elastomeric Threads.....	202
Archwire Protection Sleeves	203



Elastomeric Color Chains

Available in a wide range of colors that remain stable with time, our color chains are sold in 4.5 meter spools in either long or short.

Color		Short	Long
	Silver	88-1001-0000	88-3001-0000
	Pearl	88-1002-0000	—
	Clear	88-1003-0000	88-3003-0000
	Blue	88-1004-0000	88-3004-0000
	Red	88-1005-0000	—
	Green	88-1006-0000	88-3006-0000
	Black	88-1007-0000	88-3007-0000
	White	88-1008-0000	—
	Pink	88-1009-0000	—
	Metallic Purple	88-1010-0000	—
	Metallic Blue	88-1011-0000	—
	Metallic Red	88-1012-0000	—
	Lavender	88-1013-0000	88-3013-0000
	Fuchsia	88-1014-0000	—
	Yellow	88-1015-0000	—
	Orange	88-1016-0000	88-3016-0000
	Navy Blue	88-1017-0000	—
	Light Blue	88-1018-0000	88-3018-0000
	Maroon	88-1019-0000	—
	Grey	88-1020-0000	88-3020-0000
	Glow Blue	88-1021-0000	—
	Metallic Gold	88-1022-0000	88-3022-0000

Orchestra Intra and Extraoral Elastics

Our range of intra and extraoral elastics are made from the highest quality medical grade latex. The stated force corresponds to the force three times its original listed diameter. Intraoral elastics are packaged by 5000 pcs. Extraoral elastics are packaged by 2500 pcs. Each bag contain 1 insertion plastic instrument. Latex free elastics are listed on the chart below with product numbers (P/N) in white.

EXTRA-ORAL				
Force Size	2.5 oz (70,88 g)	4.5 oz (113,40 g)	6.5 oz (170,10 g)	8,5 oz (240, 98 g)
3 mm (1/8")	85-0321-0000 Harp	85-0331-0000 Xylophone	85-0341-0000 Tuba	
4 mm (3/16")	85-0421-0000 Tambourine	85-0431-0000 Piano	85-0441-0000 Trumpet	85-0461-0000 Trombone
4 mm (3/16") (Latex Free)		85-0430-0000 Jazz (Latex Free)	85-0440-0000 Rock-n-Roll (Latex Free)	
6 mm (1/4")	85-0621-0000 Oboe	85-0631-0000 Flute	85-0641-0000 Baritone	85-0661-0000 French Horn
6 mm (1/4") (Latex Free)		85-0630-0000 Hip-Hop (Latex Free)	85-0640-0000 Classical (Latex Free)	
8 mm (5/16")	85-0821-0000 Cello	85-0831-0000 Clarinet	85-0841-0000 Saxophone	85-0861-0000 Guitar
10 mm (3/8")	85-1021-0000 Cymbals	85-1031-0000 Tympani		85-1061-0000 Violin



GC Ties Elastomeric Ligatures

Latex-free, GC Ties ligatures allow for an upper and lower archwire change with a single stick. The high quality of the polyurethane enables these ligatures to recover 90 % of their initial shape after being stretched.

Available in bags of 50 sticks, 20 ligatures per stick. Outer diameter: .120"

■	Silver	87-1201-0000	■	Lavender	87-1213-0000
□	Pearl	87-1202-0000	■	Fuchsia	87-1214-0000
□	Clear	87-1203-0000	■	Yellow	87-1215-0000
■	Blue	87-1204-0000	■	Orange	87-1216-0000
■	Red	87-1205-0000	■	Navy Blue	87-1217-0000
■	Green	87-1206-0000	■	Light Blue	87-1218-0000
■	Black	87-1207-0000	■	Maroon	87-1219-0000
□	White	87-1208-0000	■	Grey	87-1220-0000
■	Pink	87-1209-0000	■	Glow Blue	87-1221-0000
■	Metallic Purple	87-1210-0000	■	Golden	87-1222-0000
■	Metallic Blue	87-1211-0000	■	Glow Clear	87-1223-0000
■	Metallic Red	87-1212-0000	■	Tooth	87-1224-0000

Separators



Very strong, latex-free separators with memory effect, these separators are ideal for creating spaces between molars. Sold by 1000.

Description	P/N
Radiopaque Separators	89-1602-0000



Hollow or Solid Elastomeric Threads

Our hollow or solid elastomeric threads are available in diameter sizes of .025" and .030" in 15 m spools.

Hollow Threads	P/N
.025" Grey	89-4325-0000
.025" Clear	89-3325-0000
.030" Clear	89-3335-0000

Solid Threads	P/N
.025" Grey	89-4125-0000
.025" Clear	89-3125-0000
.030" Grey	89-4135-0000
.030" Clear	89-3135-0000

Solid Fiber with Treated Surface	P/N
.025" Clear	89-3222-0000
.025" Grey	89-4222-0000
.030" Grey	89-4232-0000



Archwire Protection Sleeves

Available in 3 m Spools

For Archwires up to .018"		P/N
Grey		89-4311-0000
Clear		89-3311-0000

For Archwires up to .021" x .027"		P/N
Grey		89-4321-0000
Clear		89-3321-0000

For .045" Archwires		P/N
Grey		89-4341-0000

