

Bands



Global Experts in Orthodontics

03

**GC Orthodontics**  
GLOBAL EXPERTS IN ORTHODONTICS

## Bands

Sure-Snap Molar Bands	85
Sure-Snap Bicuspid Bands	86
A-Fit Molar Bands	87
Legacy Molar Bands	88
Weldable Lingual Attachments	89
Band Boxes	90





# Sure-Snap Molar Bands

Sure-Snap Molar Bands are made from a special alloy creating a band that offers at the same time solidity and elasticity. The difference in thickness between occlusal and gingival edges creates bands that are:

- Resistant to the pressure of the band-pusher (occlusal edge)
- Elastic enough to fit over the contour of the tooth without deformation (gingival edge)

Sure-Snap features a unique anatomical form for each type of molar, an easy-to-read permanent identification and an electroeroded interior for increased bond strength.

These easy-to-fit bands are available in 32 sizes for first molars and 30 sizes for second molars.



## Sure-Snap First Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2									
Intro	1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	3	3	3	3	2	2	2	2	2	1	1	1	1	1	1	1	
Standard	1	1	1	1	2	2	2	3	4	4	4	5	5	6	6	6	6	6	6	5	5	4	3	2	2	2	1	1	1	1	1	
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	14	12	10	8	6	4	4	2	2	2	2	1	1	1

### Upper First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 32 Sizes  
Standard Kit: 200 Assorted Bands, 32 Sizes  
Master Kit: 400 Assorted Bands, 32 Sizes



UL18

**Dental arch**

U Upper (maxillary)  
L Lower (mandibular)

**Side**

L Left  
R Right

**Size**

Sizes range from 1 to 32

### Lower First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 32 Sizes  
Standard Kit: 200 Assorted Bands, 32 Sizes  
Master Kit: 400 Assorted Bands, 32 Sizes

## Sure-Snap Second Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2									
Intro	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	3	3	3	3	2	2	2	2	2	2	1	1	1	1	1		
Standard	1	1	1	2	2	2	2	3	4	4	4	5	5	6	6	6	6	6	6	5	5	4	3	2	2	2	2	1	1	1		
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	16	12	10	8	6	4	4	2	2	2	2	1		

### Upper Second Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 30 Sizes  
Standard Kit: 200 Assorted Bands, 30 Sizes  
Master Kit: 400 Assorted Bands, 30 Sizes



SUL26

**Molar**

S Means second molar

**Dental arch**

U Upper (maxillary)  
L Lower (mandibular)

**Side**

L Left  
R Right

**Size**

Sizes range from 1 to 30

### Lower Second Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 30 Sizes  
Standard Kit: 200 Assorted Bands, 30 Sizes  
Master Kit: 400 Assorted Bands, 30 Sizes

# Sure-Snap Bicuspid Bands


The same high quality manufacturing that goes into our molar bands is in our bicuspid bands as well. Band seating is accomplished through a rocking movement, alternating mesial and distal pressure. Available in three anatomical forms:

- Maxillary First Bicuspid (32 sizes)
- Mandibular First Bicuspid (25 sizes)
- Mandibular Second Bicuspid (27 sizes)

Thanks to their easy seating, you will not hesitate to use Sure-Snap Bicuspid Bands instead of brackets for severe rotations or in case of bond failure.

## Upper Bicuspid, Bands Only

Starter Kit: 50 Assorted Bands, 22 Sizes  
Intro Kit: 100 Assorted Bands, 32 Sizes  
Standard Kit: 200 Assorted Bands, 32 Sizes  
Master Kit: 400 Assorted Bands, 32 Sizes




### Sure-Snap Upper Bicuspid Kits – Distribution

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
Starter						2	2	2	2	2	2	2	2	2	2	4	4	4	2	2	2	2	2	2	2	2	2						
Intro	2	2	2	2	2	2	2	2	2	2	4	4	4	4	6	6	6	6	6	4	4	4	4	2	2	2	2	2	2	2	2	2	2
Standard	2	2	2	2	2	4	6	6	6	6	6	8	8	10	12	14	14	14	12	10	8	8	6	6	6	4	4	4	2	2	2	2	
Master	4	4	4	6	6	10	10	14	14	16	16	16	18	18	20	22	22	22	20	18	16	16	14	14	12	12	8	8	8	4	4	4	

## 1st Lower Bicuspid, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 25 Sizes  
Standard Kit: 200 Assorted Bands, 25 Sizes  
Master Kit: 400 Assorted Bands, 25 Sizes



### Sure-Snap First Lower Bicuspid Kits - Distribution

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
Starter			2	4	4	4	4	4	4	4	4	4	4	4	4										
Intro	2	2	2	2	2	2	4	4	4	4	4	6	6	8	8	8	8	4	4	4	4	2	2	2	2
Standard	2	2	2	2	4	4	8	8	8	8	10	12	12	18	18	18	18	10	8	8	8	4	4	2	2
Master	4	4	4	6	6	8	8	14	14	18	18	24	24	36	36	36	36	24	24	16	16	8	8	4	4

## 2nd Lower Bicuspid, Bands Only

Starter Kit: 50 Assorted Bands, 27 Sizes  
Intro Kit: 100 Assorted Bands, 27 Sizes  
Standard Kit: 200 Assorted Bands, 27 Sizes  
Master Kit: 400 Assorted Bands, 27 Sizes



### Sure-Snap Second Lower Bicuspid Kits - Distribution

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter						1	1	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	1	1
Intro						2	2	2	2	2	2	2	4	4	4	6	6	8	8	8	6	6	4	4	4	2	2	2	2	2	2	2
Standard						2	2	2	4	4	6	6	6	8	10	12	14	16	16	16	14	12	10	8	6	6	6	4	4	2	2	2
Master						2	4	4	6	8	10	12	12	16	20	26	28	34	34	34	28	26	20	16	14	12	10	8	6	4	4	2



# A-Fit Molar Bands

A-Fit bands are more formable, medium tempered bands. With a slightly higher crown height, they have the strength to withstand seating and functional pressure, but are adaptive to the tooth's contour as well. A-Fit bands are marked for easy size identification and etched for increased bond strength. Available in an even broader range of sizes for your added convenience.



## A-Fit First Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)

Size	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	
Starter									1	2	2	2	2	2	2	2	2	2	2	2	2														
Intro	1	1	1	1	1	1	1	1	1	2	2	2	2	2	3	3	3	2	2	2	2	2	1	1	1	1	1	1	1	1	1	1	1	1	
Std.	1	1	1	1	2	2	3	4	4	4	5	6	6	6	6	6	5	5	4	4	3	3	2	1	1	1	1	1	1	1	1	1	1	1	
Master	2	2	2	2	4	4	6	8	8	8	10	12	12	12	12	12	10	10	8	8	6	6	4	2	2	2	2	2	2	2	2	2	2	2	

### Upper First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 34 Sizes  
Standard Kit: 200 Assorted Bands, 34 Sizes  
Master Kit: 400 Assorted Bands, 34 Sizes



25UL

### Lower First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 34 Sizes  
Standard Kit: 200 Assorted Bands, 34 Sizes  
Master Kit: 400 Assorted Bands, 34 Sizes

**Size**  
Sizes range from 5 to 38

**Dental arch**  
U Upper (maxillary)  
L Lower (mandibular)

**Side**  
L Left  
R Right

## A-Fit Second Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
Starter											1	2	2	2	2	2	2	2	2	2	2	2											
Intro	1	1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	3	3	3	3	2	2	2	2	1	1	1	1	1	1	1	1
Standard	1	1	1	1	2	2	2	3	4	4	4	5	5	6	6	6	6	6	5	5	4	3	2	2	2	1	1	1	1	1	1	1	1
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	14	12	10	8	6	4	4	2	2	2	2	1	1	1	1

### Upper Second Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 32 Sizes  
Standard Kit: 200 Assorted Bands, 32 Sizes  
Master Kit: 400 Assorted Bands, 32 Sizes



1SUL

### Lower Second Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
Intro Kit: 100 Assorted Bands, 32 Sizes  
Standard Kit: 200 Assorted Bands, 32 Sizes  
Master Kit: 400 Assorted Bands, 32 Sizes

**Size**  
Sizes range from 1 to 32

**Molar**  
S Means second molar

**Dental arch**  
U Upper (maxillary)  
L Lower (mandibular)

**Side**  
L Left  
R Right

# Legacy Molar Bands

Legacy first molar bands were designed with a shorter crown height, per the request of a number of orthodontists. The temper is softer, making sizing, placement and adaptation simple. Legacy bands hold up well throughout treatment. They are available in 32 sizes.



## Legacy First Molar Kits – Distribution by Quadrant (UR,UL,LR,LL)


Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Starter											1	2	2	2	2	2	2	2	2	2	2	2	2									
Intro	1	1	1	1	1	1	1	1	1	1	2	2	2	2	2	2	2	3	3	3	3	2	2	2	1	1	1	1	1	1	1	1
Standard	1	1	1	1	2	2	2	3	4	4	4	5	5	6	6	6	6	6	6	5	5	4	3	2	2	2	1	1	1	1	1	1
Master	1	1	1	2	2	2	4	4	6	6	8	10	10	12	14	16	16	16	14	12	10	8	6	4	4	2	2	2	2	1	1	1

### Upper First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
 Intro Kit: 100 Assorted Bands, 32 Sizes  
 Standard Kit: 200 Assorted Bands, 32 Sizes  
 Master Kit: 400 Assorted Bands, 32 Sizes

### Lower First Molars, Bands Only

Starter Kit: 50 Assorted Bands, 13 Sizes  
 Intro Kit: 100 Assorted Bands, 32 Sizes  
 Standard Kit: 200 Assorted Bands, 32 Sizes  
 Master Kit: 400 Assorted Bands, 32 Sizes



## GC Orthodontics Band Configurations

At GC Orthodontics, we are proud of our reputation for high quality products designed to meet the needs of an ever-changing orthodontic world. We have a variety of band, tube and lingual attachment combinations that are manufactured and assembled with the same passion for excellence as our other premium products. If you don't see it, don't assume we can not provide it. Ask your GC Orthodontics representative for assistance. We are proud to be Global Experts in Orthodontics!



Upper 1st Molar  
Triple Convertible  
Lingual Cleat



Upper 1st Molar  
Double Convertible  
Lingual Sheath



Upper 1st Molar  
Single Convertible  
Lingual Sheath



Lower 1st Molar  
Double Convertible  
Lingual Cleat




Lower 1st Molar  
Single Convertible  
Lingual Seating Lug


# Weldable Lingual Attachments




## Lingual Attachments

		PART NUMBER	
		RIGHT	LEFT
	Standard, with locking indent	43-0250-0000	
	Cetlin 8° offset with V-Stop	43-027R-0000	43-027L-0000
	With hook and locking indent	43-028R-0000	43-028L-0000
	Without locking indent	43-0260-0000	

		PART NUMBER
	Weldable Button, Small Flat	43-0310-0000
	Weldable Button, Medium Curved	43-0300-0000
	Weldable Button, Medium Flat	43-0320-0000

		PART NUMBER
	Weldable Seating Lug, Flat	43-0340-0000
	Weldable Seating Lug, Curved	43-0350-0000

		PART NUMBER
	Weldable Cleat, Small (3.3mm)	43-0210-0000
	Weldable Cleat, Medium (4.4mm)	43-0220-0000
	Weldable Cleat, Large (5.2mm)	43-0230-0000
	Weldable Cleat, Extra Large (8.4mm)	43-0240-0000

		PART NUMBER
	Weldable Ball Hook	43-0380-0000

# Band Boxes

GC Orthodontics' band boxes are designed to help you keep your bands safe and organized, as well as to enhance the decor of an orthodontic practice. Our band boxes are color coordinated for easy identification, and they travel well.

