


## Bondable Tubes and Lingual Attachments



Global Experts in Orthodontics

03

**GC Orthodontics**  
GLOBAL EXPERTS IN ORTHODONTICS

## Bondable Tubes and Lingual Attachments

---

Experience Self-Ligating Upper & Lower 1st Molar Tubes	61
Legend Low Profile 1st & 2nd Molar Tubes	62
Nonconvertible Upper & Lower 1st Molar Tubes	64
Convertible Upper 1st Molar Tubes	65
Convertible Lower 1st Molar Tubes	66
Upper 2nd Molar Tubes	67
Lower 2nd Molar Tubes	68
Access 2nd Molar Tubes and Mini Eruption Tubes	69
Lingual Attachments	70

---

# EXPERIENCE




## EXPERIENCE self-ligating buccal tubes



This tube complements the EXPERIENCE appliances on 1st molars. It features all the benefits of EXPERIENCE brackets:




- Active clip easy to open from the buccal face and easy to close at the fingertip.
- Very large micro-sandblasted base that covers the entire buccal surface of the molar.
- Low profile for patient comfort and less bond failure.
- Available in several prescriptions
- MBT\* is available with standard or small base (SB)




### Standard Base


#### Maxillary Single Self-Ligating 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.0	-10°	14°D	UR6	20-1160-0010	20-2160-0010	Roth
				UL6	20-1260-0010	20-2260-0010	
	4.0	-14°	10°D	UR6	20-1160-0020	20-2160-0020	MBT
				UL6	20-1260-0020	20-2260-0020	

#### Mandibular Single Self-Ligating 1st Molar


	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.0	-25°	4°D	LR6	20-1460-0010	20-2460-0010	Roth
				LL6	20-1360-0010	20-2360-0010	
	4.0	-20°	0°D	LR6	20-1460-0020	20-2460-0020	MBT
				LL6	20-1360-0020	20-2360-0020	
	4.0	-10°	4°D	LR6	—	20-2460-0032	XT
				LL6	—	20-2360-0032	

#### Maxillary and Mandibular Double Self-Ligating 1st Molar


	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.0	-20°	20°D	UR6	—	20-2160-0030	XT
				UL6	—	20-2260-0030	
	4.0	-10°	4°D	LR6	—	20-2460-0030	XT
				LL6	—	20-2360-0030	

### Small Base

#### Maxillary Single Self-Ligating 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.0	-14°	10°D	UR6	20-1160-S020	20-2160-S020	MBT
				UL6	20-1260-S020	20-2260-S020	

#### Mandibular Single Self-Ligating 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.0	-20°	0°D	LR6	20-1460-S020	20-2460-S020	MBT
				LL6	20-1360-S020	20-2360-S020	

# LEGEND LP

## LEGEND®

### LOW PROFILE BUCCAL TUBES

LEGEND'S LOW PROFILE HEIGHT AND MALLEABLE HOOK DESIGN BOOSTS PATIENT COMFORT, WHILE EASY TO GRIP POSITIONING GUIDES AID IN THE SECURE HANDLING AND EASY PLACEMENT IN THE POSTERIOR DURING BONDING. COUPLED WITH MICRO-ETCHED 80-GAUGE MESH AND ANATOMICAL CONTOURED BONDING PADS, LEGEND BUCCAL TUBES DELIVER MAXIMUM BOND STRENGTH AND SUPERIOR PERFORMANCE




#### Legend Low Profile Non-Convertible Bondable Buccal Tubes - 1st Molar


Color Codes



##### Maxillary Single Low Profile Non-Convertible 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
						4.5	
				UL6	40-1260-0114	40-2260-0114	
	4.5	-14°	10°D	UR6	40-1160-0106	40-2160-0106	MBT
				UL6	40-1260-0106	40-2260-0106	
	4.5	-14°	0°D	UR6	40-1160-0107	40-2160-0107	Roth/MBT (Opt)
				UL6	40-1260-0107	40-2260-0107	
	4.5	-10°	13°D	UR6	40-1160-0115	40-2160-0115	Alex
				UL6	40-1260-0115	40-2260-0115	
	4.5	-10°	8°D	UR6	40-1160-0116	40-2160-0116	And
				UL6	40-1260-0116	40-2260-0116	

##### Mandibular Single Low Profile Non-Convertible 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
						4.5	
				LL6	40-1360-0114	40-2360-0114	
	4.5	-20°	0°D	LR6	40-1460-0106	40-2460-0106	MBT
				LL6	40-1360-0106	40-2360-0106	
	4.5	-10°	8°D	LR6	40-1460-0115	40-2460-0115	Alex
				LL6	40-1360-0115	40-2360-0115	
	4.5	-25°	0°D	LR6	40-1460-0116	40-2460-0116	And
				LL6	40-1360-0116	40-2360-0116	



# LEGEND<sup>®</sup>

## LOW PROFILE BUCCAL TUBES




### Legend Low Profile Bondable Buccal Tubes - 2nd Molar


Color Codes



#### Maxillary Single Low Profile 2nd Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
						3.6	
				UL7	40-1270-0114	40-2270-0114	
	3.6	-14°	10°D	UR7	40-1170-0106	40-2170-0106	MBT
				UL7	40-1270-0106	40-2270-0106	
	3.6	-14°	0°D	UR7	40-1170-0107	40-2170-0107	Roth/MBT (Opt)
				UL7	40-1270-0107	40-2270-0107	
	3.6	-10°	10°D	UR7	40-1170-0115	40-2170-0115	Roth/Alex
				UL7	40-1270-0115	40-2270-0115	
	3.6	-20°	10°D	UR7	40-1170-0120	40-2170-0120	GCO
				UL7	40-1270-0120	40-2270-0120	
	3.6	-10°	7°D	UR7	40-1170-0116	40-2170-0116	And
				UL7	40-1270-0116	40-2270-0116	

#### Mandibular Single Low Profile 2nd Molar


	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
						3.6	
				LL7	40-1370-0114	40-2370-0114	
	3.6	-10°	0°D	LR7	40-1470-0106	40-2470-0106	MBT
				LL7	40-1370-0106	40-2370-0106	
	3.6	0°	5°D	LR7	40-1470-0115	40-2470-0115	Alex
				LL7	40-1370-0115	40-2370-0115	
	3.6	-20°	0°D	LR7	40-1470-0120	40-2470-0120	GCO
				LL7	40-1370-0120	40-2370-0120	
	3.6	-25°	0°D	LR7	40-1470-0116	40-2470-0116	And
				LL7	40-1370-0116	40-2370-0116	

# Non-Convertible Tubes


## Bondable Non-Convertible Buccal Tubes - 1st Molar

### VLP Low Profile Buccal Tubes


#### Maxillary Single Non-Convertible 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.2	-14°	10°D	UR6	40-1160-0120	40-2160-0120	Roth
			UL6	40-1260-0120	40-2260-0120		
4.2	-10°	8°D	UR6	40-1160-0108	40-2160-0108	MBT	
			UL6	40-1260-0108	40-2260-0108		

#### Mandibular Single Non-Convertible 1st Molar


	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.2	-25°	5°D	LR6	40-1460-0105	40-2460-0105	Roth
			LL6	40-1360-0105	40-2360-0105		
4.2	-20°	0°D	LR6	40-1460-0120	40-2460-0120	MBT	
			LL6	40-1360-0120	40-2360-0120		

#### Maxillary Single Non-Convertible 1st Molar - Wide Entry

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	5.5	-10°	7°D	UR6	—	40-2160-0035	Andrews
			UL6	—	40-2260-0035		

## Double Non-Convertible Buccal Tubes - 1st Molar


### Mandibular Double Non-Convertible 1st Molar (Auxillary Tubes)

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
	4.0	0°	0°D	LR/LL6	40-1660-0044	40-2660-0044	Double Non-Convertible Universal Archwire and Aux. Slot Same Size	Bio


## AXCESS Tubes - Bondable

These entry point buccal tubes are manufactured in one-piece using MIM (Metal Injection Molded) technology

#### Maxillary Single Non-Convertible 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.4	-10°	6°D	UR6	40-1160-0025	40-2160-0025	Roth/And
			UL6	40-1260-0025	40-2260-0025		
4.4	0°	6°D	UR6	40-1160-0026	40-2160-0026	Edgewise	
			UL6	40-1260-0026	40-2260-0026		

#### Mandibular Single Non-Convertible 1st Molar



	Length mm	Torque	Rotation	Tooth	BONDABLE		Rx
					.018 SLOT	.022 SLOT	
	4.4	-25°	6°D	LR6	40-1460-0025	40-2460-0025	Roth/And
			LL6	40-1360-0025	40-2360-0025		
4.4	0°	6°D	LR6	40-1460-0026	40-2460-0026	Edgewise	
			LL6	40-1360-0026	40-2360-0026		



# Convertible Tubes


## Single Convertible Buccal Tubes - Upper 1st Molar

### Maxillary Single Convertible 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
	3.8	-10°	14°D	UR6	40-1160-0012	40-2160-0012	Single Convertible	Roth
				UL6	40-1260-0012	40-2260-0012		
	4.0	-14°	14°D	UR6	40-1160-0061	40-2160-0061	Single Convertible	Roth
				UL6	40-1260-0061	40-2260-0061		
	3.8	-10°	7°D	UR6	40-1160-0034	40-2160-0034	Single Convertible	And
				UL6	40-1260-0034	40-2260-0034		
3.8	-14°	10°D	UR6	40-1160-0022	40-2160-0022	Single Convertible	MBT	
			UL6	40-1260-0022	40-2260-0022			
4.2	-14°	10°D	UR6	40-1160-0110	40-2160-0110	Single Convertible	MBT	
			UL6	40-1260-0110	40-2260-0110			
	4.0	-10°	14°D	UR6	40-1160-0016	40-2160-0016	Single Convertible - Large Entry	Roth
				UL6	40-1260-0016	40-2260-0016		


## Double Convertible Buccal Tubes (No HG tube) - Upper 1st Molar

### Maxillary Single Convertible 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
	4.0	-14°	14°D	UR6	40-1160-0011	40-2160-0011	Double Convertible No HG .018 x .025 Auxillary Slot	Roth
				UL6	40-1260-0011	40-2260-0011		
	4.0	-14°	14°D	UR6	40-1160-0055	40-2160-0055	Double Convertible No HG .022 x .028 Auxillary Slot	Roth
				UL6	40-1260-0055	40-2260-0055		
	4.0	-20°	20°D	UR6	40-1160-0030	40-2160-0030	Double Convertible No HG .018 x .025 Auxillary Slot	XT
				UL6	40-1260-0030	40-2260-0030		
4.0	-20°	20°D	UR6	—	40-2160-0031	Double Conv. No HG, Large Base .018 x .025 Auxillary Slot	XT	
			UL6	—	40-2260-0031			

## Double Convertible Buccal Tubes with HG tube - Upper 1st Molar

### Maxillary Single Convertible 1st Molar








	Length mm	Torque	Rotation		Tooth	BONDABLE		Rx	Rx
			(Slot)	(HG)		.018 SLOT	.022 SLOT		
	4.0	-14°	10°D	10°D	UR6	40-1160-0020	40-2160-0020	Double Convertible .045 HG Gingival	MBT
					UL6	40-1260-0020	40-2260-0020		
	4.0	-10°	7°D	7°D	UR6	40-1160-0079	40-2160-0079	Double Convertible .045 HG Occlusal	And
					UL6	40-1260-0079	40-2260-0079		



# Convertible Tubes






## Single Convertible Buccal Tubes - Lower 1st Molar

### Mandibular Single Convertible 1st Molar

	Length mm	Torque	Ang	Rot	Tooth	BONDABLE		Details	Rx
						.018 SLOT	.022 SLOT		
	3.8	-25°	0°	4°D	LR6	40-1460-0012	40-2460-0012	Single Convertible	Roth
					LL6	40-1360-0012	40-2360-0012		
	4.0	-30°	-1°	4°D	LR6	40-1460-0062	40-2460-0062	Single Convertible	Roth
					LL6	40-1360-0062	40-2360-0062		
	4.0	-20°	0°	0°D	LR6	40-1460-0022	40-2460-0022	Single Convertible	MBT
					LL6	40-1360-0022	40-2360-0022		
	3.8	-25°	0°	0°D	LR6	40-1460-0034	40-2460-0034	Single Convertible	And
					LL6	40-1360-0034	40-2360-0034		
	3.6	-22°	0°	4°D	LR6	40-1460-0029	—	Single Convertible	Polar
					LL6	40-1360-0029	—		
	3.8	0°	0°	6°D	LR6	—	40-2460-0092	Single Convertible	Edge
					LL6	—	40-2360-0092		
	3.8	0°	0°	0°D	0°	40-1460-0093	40-2460-0093	Single Convertible	Edge
						40-1360-0093	40-2360-0093		
	3.8	-10°	0°	0°D	LR6	40-1460-0076	40-2460-0076	Single Convertible	XT
					LL6	40-1360-0076	40-2360-0076		
	4.0	-10°	0°	0°D	LR6	40-1460-0016	40-2460-0016	Single Convertible - Large Entry	XT
					LL6	40-1360-0016	40-2360-0016		

## Double Convertible Buccal Tubes - Lower 1st Molar

### Mandibular Double Convertible 1st Molar (Auxillary Tubes)

	Length mm	Torque	Ang	Rot	Tooth	BONDABLE		Details	Rx
						.018 SLOT	.022 SLOT		
	3.8	-25°	0°	4°D	LR6	40-1460-0011	40-2460-0011	Double Convertible .018 x .025 Auxillary Slot	Roth
					LL6	40-1360-0011	40-2360-0011		
	4.0	-30°	-1°	4°D	LR6	—	40-2460-0060	Double Convertible .022 x .028 Auxillary Slot	Roth
					LL6	—	40-2360-0060		
	3.6	-10°	0°	4°D	LR6	—	40-2460-0031	Double Convertible Wide Base .018 x .025 Auxillary Slot	XT
					LL6	—	40-2360-0031		
	4.0	-30°	-1°	4°D	LR6	—	40-2460-0058	Double Convertible .022 x .028 Aux Slot Occlusal Hook	Roth
					LL6	—	40-2360-0058		
	4.0	0°	0°	0°D	LR6	40-1660-0045	40-2660-0045	Double Convertible Universal .018 x .025 Auxillary Slot	Bio
					LL6				





# 2nd Molar Tubes

## Single Non-Convertible Buccal Tubes - Upper 2nd Molar

### Maxillary 2nd Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
	4.0	-14°	14°D	UR7	40-1170-0013	40-2170-0013	Single Non-Convertible Distal Extension	Roth
				UL6	40-1270-0013	40-2270-0013		
	3.0	-10°	14°D	UR7	—	40-2170-0039	Single Non-Convertible, Small Distal Extension	Roth
				UL7	—	40-2270-0039		
	3.0	-10°	7°D	UR7	40-1170-0048	40-2170-0048	Single Non-Convertible, Small	And
				UL7	40-1270-0048	40-2270-0048		
	3.0	-10°	7°D	UR7	40-1170-0049	40-2170-0049	Single Non-Convertible, Small Wide Base	And
				UL7	40-1270-0049	40-2270-0049		
	4.0	-14°	14°D	UR7	40-1170-0014	40-2170-0014	Single Non-Convertible	Roth
				UL7	40-1270-0014	40-2270-0014		
	4.0	-14°	0°D	UR7	40-1170-0063	40-2170-0063	Single Non-Convertible	Roth
				UL7	40-1270-0063	40-2270-0063		
	4.0	-20°	20°D	UR7	40-1170-0032	40-2170-0032	Single Non-Convertible	XT
				UL7	40-1270-0032	40-2270-0032		
	4.0	-14°	10°D	UR7	40-1170-0022	40-2170-0022	Single Non-Convertible	MBT
				UL7	40-1270-0022	40-2270-0022		
	4.0	-10°	7°D	UR7	40-1170-0034	40-2170-0034	Single Non-Convertible	And
				UL7	40-1270-0034	40-2270-0034		
	3.8	-14°	14°D	UR7	40-1170-0012	40-2170-0012	Single Non-Convertible	Roth
				UL7	40-1270-0012	40-2270-0012		
	3.5	-10°	7°D	UR7	40-1170-0018	40-2170-0018	Single Non-Convertible Peucat Wide Entry	And
				UL7	40-1270-0018	40-2270-0018		

# 2nd Molar Tubes

## Single Non-Convertible Buccal Tubes - Lower 2nd Molar

### Mandibular 2nd Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
	3.0	-25°	4°D	LR7	40-1470-0049	40-2470-0049	Single Non-Convertible, Small Distal Extension, Larger Pad	Roth
				LL7	40-1370-0049	40-2370-0049		
	3.0	-25°	4°D	LR7	40-1470-0048	40-2470-0048	Single Non-Convertible, Small Distal Extension	Roth
				LL7	40-1370-0048	40-2370-0048		
	4.0	-25°	6°D	LR7	40-1470-0012	40-2470-0012	Single Non-Convertible	Roth
				LL7	40-1370-0012	40-2370-0012		
	4.0	-30°	4°D	LR7	40-1470-0064	40-2470-0064	Single Non-Convertible -1 Angulation	Roth
				LL7	40-1370-0064	40-2370-0064		
	4.0	-25°	0°D	LR7	40-1470-0034	40-2470-0034	Single Non-Convertible	And
				LL7	40-1370-0034	40-2370-0034		
	4.0	0°	0°D	LR7	40-1470-0095	40-2470-0095	Single Non-Convertible Distal Extension	Edge
				LL7	40-1370-0095	40-2370-0095		
	4.0	0°	6°D	LR7	40-1470-0096	40-2470-0096	Single Non-Convertible Distal Extension	Edge
				LL7	40-1370-0096	40-2370-0096		
	3.2	-10°	0°D	LR7	40-1470-0022	40-2470-0022	Single Non-Convertible	MBT
				LL7	40-1370-0022	40-2370-0022		
	3.2	0°	6°D	LR7	40-1470-0097	40-2470-0097	Single Non-Convertible	Edge
				LL7	40-1370-0097	40-2370-0097		
	3.2	0°	0°D	LR7	40-1470-0098	40-2470-0098	Single Non-Convertible	Edge
				LL7	40-1370-0098	40-2370-0098		
	4.6	0°	0°D	LR7	40-1470-0094	40-2470-0094	Single Non-Convertible	Edge
				LL7	40-1370-0094	40-2370-0094		
	3.5	-25°	4°D	LR7	40-1470-0018	40-2470-0018	Single Non-Convertible Peucat Wide Entry	Roth
				LL7	40-1370-0018	40-2370-0018		




# 2nd Molar Tubes

## AXCESS Tubes - Bondable 2nd Molar



### Maxillary 2nd Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
					4.4	-10°		
UL7	40-1270-0025	40-2270-0025						
4.0	-10°	6°D	UR7	40-1170-0088	40-2170-0088	Access Single Non-Convertible	Roth And	
			UL7	40-1270-0088	40-2270-0088			
4.0	-10°	0°D	UR7	40-1170-0091	40-2170-0091	Access Single Non-Convertible	Roth And	
			UL7	40-1270-0091	40-2270-0091			
4.0	0°	0°D	UR7	40-1170-0090	40-2170-0090	Access Single Non-Convertible	Edge	
			UL7	40-1270-0090	40-2270-0090			

### Mandibular 2nd Molar


	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
					4.4	-25°		
LL7	40-1370-0025	40-2370-0025						
4.0	-25°	6°D	LR7	—	40-2470-0089	Access Single Non-Convertible	Roth And	
			LL7	—	40-2370-0089			
4.0	-25°	0°D	LR7	40-1470-0091	40-2470-0091	Access Single Non-Convertible	Roth And	
			LL7	40-1370-0091	40-2370-0091			
4.0	0°	0°D	LR7	40-1470-0090	40-2470-0090	Access Single Non-Convertible	Edge	
			LL7	40-1370-0090	40-2370-0090			

## Mini Eruption Tubes/Pix Tubes


(2mm offset buccal tubes are designed to fit the mesio-buccal cusp of partially erupted 2nd molars)




### Mini Eruption Tubes - Maxillary 2nd Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
					2.0	-10°		
UL7	40-1270-0019	40-2270-0019						
2.0	-10°	0°D	UR7	40-1170-0112	40-2170-0112	Mini Eruption Tube/No Hook	Roth And	
			UL7	40-1270-0112	40-2270-0112			

### Mini Eruption Tubes - Mandibular 2nd Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
					2.0	-25°		
LL7	40-1370-0019	40-2370-0019						
2.0	-25°	0°D	LR7	40-1470-0112	40-2470-0112	Mini Eruption Tube/No Hook	Roth And	
			LL7	40-1370-0112	40-2370-0112			

### Pix Tubes - Maxillary/Mandibular 2nd Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE		Details	Rx
					.018 SLOT	.022 SLOT		
					2.0	0°		
UL7								

# Bondable Lingual Attachments



## Lingual Attachments

Bondable Lingual Attachments	Quantity	PART NUMBER
------------------------------	----------	-------------

### Bondable Buttons

	Bondable Button, Curved Base	10/pkg	42-0010-0000
	Bondable Button, Flat Base	10/pkg	42-0020-0000
	Bondable Ceramic Button, Flat Base	10/pkg	44-1010-0001
	Bondable Seagull Button, Open	10/pkg	44-2130-0000
	Bondable Seagull Button, Closed	10/pkg	44-2140-0000



### Bondable Cleats

	Bondable Cleat, Curved Base	10/pkg	42-0050-0000
	Bondable Cleat, Flat Base	10/pkg	42-0060-0000

### Bondable Eyelets

	Bondable Eyelet	10/pkg	42-0070-0000
---	-----------------	--------	--------------


### Bondable Lingual Loops/Surgical Buttons

	Lingual Loop For Curved Teeth	10/pkg	42-0030-0000
	Lingual Loop For Flat Teeth	10/pkg	42-0040-0000

## Comfortube Universal 2nd Molar

Designed for lingual but also useful for some labial applications

### Maxillary Single Non-Convertible 1st Molar

	Length mm	Torque	Rotation	Tooth	BONDABLE	Details	Rx
	5.0	0°	0°D	R/L	.018 x .025 SLOT		
					40-1770-0000	Single Non-Convertible	Lingual