



Elastics and Elastomerics



Global Experts in Orthodontics

03

GC Orthodontics
GLOBAL EXPERTS IN ORTHODONTICS



















Elastics and Elastomerics

Elastics	113
Elastomerics	114
Colored Elastomeric Ligature Ties	114
Colored Elastomeric Chain	115
GC Ortho Super Chain	115
Separators	115
Elastomeric Thread	115
Archwire Sleeve/Bumper Tubing	115
Rotation Wedges	115

ORCHESTRA ELASTICS









LATEX ELASTICS

SIZE \ FORCE		LIGHT	MEDIUM	HEAVY	EXTRA HEAVY
		2.5 oz	4.0 oz	6.0 oz	8.5 oz
1/8"	LATEX	Harp 	Xylophone 	Tuba 	
		PART NUMBER	85-0321-0000	85-0331-0000	85-0341-0000
3/16"	LATEX	Tambourine 	Piano 	Trumpet 	Trombone 
		PART NUMBER	85-0421-0000	85-0431-0000	85-0441-0000
1/4"	LATEX	Oboe 	Flute 	Baritone 	French Horn 
		PART NUMBER	85-0621-0000	85-0631-0000	85-0641-0000
5/16"	LATEX	Cello 	Clarinet 	Saxophone 	Guitar 
		PART NUMBER	85-0821-0000	85-0831-0000	85-0841-0000
3/8"	LATEX	Cymbals 	Tympani 		Violin 
		PART NUMBER	85-1021-0000	85-1031-0000	

*All latex elastics are 5,000 total elastics per package (50 patient packs with 100 elastics per pack)







NEON COLORED ELASTICS

SIZE \ FORCE		MEDIUM	HEAVY
		4 oz	6 oz
3/16"	NEON COLORED LATEX	Piano 	Trumpet 
		PART NUMBER	85-0431-1111
1/4"	NEON COLORED LATEX	Flute 	Baritone 
		PART NUMBER	85-0631-1111
5/16"	NEON COLORED LATEX	Clarinet 	Saxophone 
		PART NUMBER	85-0831-1111

*All Neon Colored and Non-Latex elastics are 5,000 total elastics per package (50 patient packs with 100 elastics per pack)



NON-LATEX ELASTICS

SIZE \ FORCE		MEDIUM	HEAVY
		4 oz	6 oz
3/16"	NON-LATEX	Jazz 	Rock-n-Roll 
		PART NUMBER	85-0430-0000
1/4"	NON-LATEX	Hip-Hop 	Classical 
		PART NUMBER	85-0630-0000
5/16"	NON-LATEX	COUNTRY 	Reggae 
		PART NUMBER	85-0830-0000

Colored Ligatures and Chain

GC Colored Ligatures and Chain



GC 10-Ties and Long Stix Ligature Ties

GC 10-Ties are packaged in packs of 100 sticks (1000 ties) and Long Stix ligature ties are in packs of 42 sticks (1,008 ties). Made from elastomeric polyurethane, Our ligatures are latex-free and contain no BPA, BPF, BPS, DEHP, or plasticizers.

GC Colored Elastomeric Chain



All GC Colored Elastomeric Chain is packaged in 15' spools. GC colored elastomeric chain has optimum strength and retention, and is available in closed, short and long filaments. Our chain is color stable and resists staining. All GC colored elastomeric chain is latex-free and contain no BPA, BPF, BPS, DEHP, or plasticizers.

GC COLORS						
		PART NUMBERS				
		GC 10-Ties	GC Long Stix	Closed Chain	Short Chain	Long Chain
SILVER	87-2201-0000	87-1201-0000	88-4001-0000	88-5001-0000	88-6001-0000	
PEARL	87-2202-0000	87-1202-0000	88-4002-0000	88-5002-0000	88-6002-0000	
CLEAR	87-2203-0000	87-1203-0000	88-4003-0000	88-5003-0000	88-6003-0000	
BLUE	87-2204-0000	87-1204-0000	88-4004-0000	88-5004-0000	88-6004-0000	
RED	87-2205-0000	87-1205-0000	88-4005-0000	88-5005-0000	88-6005-0000	
GREEN	87-2206-0000	87-1206-0000	88-4006-0000	88-5006-0000	88-6006-0000	
BLACK	87-2207-0000	87-1207-0000	88-4007-0000	88-5007-0000	88-6007-0000	
WHITE	87-2208-0000	87-1208-0000	88-4008-0000	88-5008-0000	88-6008-0000	
BABY PINK	87-2209-0000	87-1209-0000	88-4009-0000	88-5009-0000	88-6009-0000	
PURPLE	87-2210-0000	87-1210-0000	88-4010-0000	88-5010-0000	88-6010-0000	
METALLIC BLUE	87-2211-0000	87-1211-0000	88-4011-0000	88-5011-0000	88-6011-0000	
METALLIC RED	87-2212-0000	87-1212-0000	88-4012-0000	88-5012-0000	88-6012-0000	
METALLIC GREEN	87-2213-0000	87-1213-0000	88-4013-0000	88-5013-0000	88-6013-0000	
FUCHSIA	87-2214-0000	87-1214-0000	88-4014-0000	88-5014-0000	88-6014-0000	
YELLOW	87-2215-0000	87-1215-0000	88-4015-0000	88-5015-0000	88-6015-0000	
ORANGE	87-2216-0000	87-1216-0000	88-4016-0000	88-5016-0000	88-6016-0000	
NAVY BLUE	87-2217-0000	87-1217-0000	88-4017-0000	88-5017-0000	88-6017-0000	
LIGHT BLUE	87-2218-0000	87-1218-0000	88-4018-0000	88-5018-0000	88-6018-0000	
MAROON	87-2219-0000	87-1219-0000	88-4019-0000	88-5019-0000	88-6019-0000	
SMOKE	87-2220-0000	87-1220-0000	88-4020-0000	88-5020-0000	88-6020-0000	
CARIBBEAN BLUE	87-2221-0000	87-1221-0000	88-4021-0000	88-5021-0000	88-6021-0000	
MARIGOLD	87-2222-0000	87-1222-0000	88-4022-0000	88-5022-0000	88-6022-0000	
GLOW CLEAR	87-2223-0000	87-1223-0000	88-4023-0000	89-5023-0000	89-6023-0000	
LILAC	87-2224-0000	87-1224-0000	88-4024-0000	88-5024-0000	88-6024-0000	
CORAL	87-2225-0000	87-1225-0000	88-4025-0000	88-5025-0000	88-6025-0000	
AQUA	87-2226-0000	87-1226-0000	88-4026-0000	88-5026-0000	88-6026-0000	
PERIWINKLE	87-2227-0000	87-1227-0000	88-4027-0000	88-5027-0000	88-6027-0000	
DARK GREEN	87-2228-0000	87-1228-0000	88-4028-0000	88-5028-0000	88-6028-0000	
LIME	87-2229-0000	87-1229-0000	88-4029-0000	88-5029-0000	88-6029-0000	
TEAL	87-2230-0000	87-1230-0000	88-4030-0000	88-5030-0000	88-6030-0000	
HOT PINK	87-2231-0000	87-1231-0000	88-4031-0000	88-5031-0000	88-6031-0000	
GOLD	87-2232-0000	87-1232-0000	88-4032-0000	88-5032-0000	88-6032-0000	
CAROLINA BLUE	87-2233-0000	87-1233-0000	88-4033-0000	88-5033-0000	88-6033-0000	
IVORY	87-2234-0000	87-1234-0000	88-4034-0000	88-5034-0000	88-6034-0000	
GREY	87-2235-0000	87-1235-0000	88-4035-0000	88-5035-0000	88-6035-0000	
QUANTITY PER PACKAGE:	1000/PKG	1008/PKG	15' SPOOL	15' SPOOL	15' SPOOL	




GC Ortho Super Chain

Made from pure high-density and high-purity materials, this chain is precision-cut and the intervals between links are strictly identical. GC Ortho Chain exerts light and continuous forces for optimal tooth movement. Available in grey or clear 2m spools.


Color	PART NUMBER		
	Closed Chain	Short Chain	Long Chain
Clear	86-1010-0000	86-1020-0000	86-1030-0000
Grey	86-2010-0000	86-2020-0000	86-2030-0000

Separators


These separators are ideal for creating spaces between molars.

	Description	Quantity	PART NUMBER
	Large Blue Radiopaque Separators	1000/pkg	89-1703-0000


Elastomeric Thread/Plastic Tie and Tubing




Solid Elastomeric Thread - Smooth (Elastic Fiber Thread)

	Description	Quantity	PART NUMBER
	.025 Clear	50' spool	89-3125-0000
.025 Grey	50' spool	89-4125-0000	
.030 Clear	50' spool	89-3135-0000	
.030 Grey	50' spool	89-4135-0000	


Solid Elastomeric Thread - Textured Non-slip (Ultra Fiber)

	Description	Quantity	PART NUMBER
	.025 Clear	25' spool	89-3222-0000
.025 Grey	25' spool	89-4222-0000	
.030 Clear	25' spool	89-3232-0000	
.030 Grey	25' spool	89-4232-0000	


Hollow Elastomeric Thread - Tie Tubing (Elast-O Fiber)

	Description	Quantity	PART NUMBER
	.025 Clear	50' spool	89-3325-0000
.025 Grey	50' spool	89-4325-0000	
.030 Clear	50' spool	89-3335-0000	
.030 Grey	50' spool	89-4335-0000	

Archwire Sleeve/Bumper Tubing (For Archwires up to .021" x .027")

	Description	Quantity	PART NUMBER
	.027 Clear	10' spool	89-3321-0000
.027 Grey	10' spool	89-4321-0000	

Rotation Wedges

	Description	Quantity	PART NUMBER
	Clear	100/pkg	69-1610-0000
Grey	100/pkg	69-1620-0000	

